



Official Website



JUMPER Health

Shenzhen Jumper Medical Equipment Co., Ltd.

Address: D Building, No. 71, Xintian Road, Fuyong Street, Baoan, Shenzhen, Guangdong, China

Tel: +86 755-26696279 **E-mail:** info@jumper-medical.com **Web:** <https://jumpermed.com>

PRODUCT BROCHURE

SINCE
1998



JUMPER

Healthcare within mind

<https://jumpermed.com>

INFRARED THERMOMETER



Forehead/Ear
Auto Switch



Magnetic
Probe Cap



Child
Mode



10s Auto
Shut Off



Smart
Fever Alert

ONE KEY MEASUREMENT, ONE SECOND READING.

Infrared Thermometer uses an infrared sensor to detect the infrared radiation emitted by the eardrum and forehead, and converts the value through a microcomputer to obtain the temperature data.

Users can quickly get measurement results after positioning properly the temperature probe in the ear or on the forehead.

In addition to adults' daily temperature measurement, it can also effectively avoid measurement difficulties caused by crying resistance when measuring infants and young children. This kind of thermometer is more suitable for the whole family to use together.

FEATURES

- Child Mode
 - Silent Mode
 - Forehead /Ear/Object Mode
- Professional Accuracy
 - Hygienic Non-Touch
 - 1 Second Reading
- 10s Auto Shut Off
 - Smart Fever Alert
 - 20 Sets Records



NEW
ARRIVAL



JPD-FR302

DISPLAY HD LCD	MEMORY STORAGE 20 SETS DATA	TEMPERATURE UNITS °C/°F
AUTOMATIC SHUTDOWN 10S ± 1S	BATTERY AAA X 2	ACCURACY(CLINICAL) 0.3°C / 0.6°F
MEASURING RANGE		
FOREHEAD MODE	22 - 43 °C / 71.6 - 109.4 °F	
EAR MODE	34 - 43 °C / 93.2 - 109.4 °F	
OBJECT MODE	0 - 100 °C / 32.0 - 212.0 °F	
Forehead/Ear Auto Switch	Magnetic Probe Cap	Child Mode
10s Auto Shut Off	Smart Fever Alert	Silent Mode

NEW
ARRIVAL



JPD-FR304

DISPLAY HD LCD	MEMORY STORAGE 20 SETS DATA	TEMPERATURE UNITS °C/°F
AUTOMATIC SHUTDOWN 10S ± 1S	BATTERY AAA X 2	ACCURACY(CLINICAL) 0.3°C / 0.6°F
MEASURING RANGE		
EAR MODE	34 - 43 °C / 93.2 - 109.4 °F	
Magnetic Probe Cap	Child Mode	10s Auto Shut Off
Smart Fever Alert	Silent Mode	10s Auto Shut Off

NEW
ARRIVAL



JPD-FR403

DISPLAY HD LCD	MEMORY STORAGE 20 SETS DATA	TEMPERATURE UNITS °C/°F
AUTOMATIC SHUTDOWN 10S ± 1S	BATTERY AAA X 2	ACCURACY(CLINICAL) 0.3°C / 0.6°F
MEASURING RANGE		
FOREHEAD MODE	22 - 43 °C / 71.6 - 109.4 °F	
EAR MODE	34 - 43 °C / 93.2 - 109.4 °F	
OBJECT MODE	0 - 100 °C / 32.0 - 212.0 °F	
Forehead/Ear Auto Switch	Magnetic Probe Cap	Child Mode
10s Auto Shut Off	Smart Fever Alert	Silent Mode

NEW ARRIVAL



JPD-FR301

DISPLAY HD LED	MEMORY STORAGE 20 SETS DATA	TEMPERATURE UNITS °C/°F
AUTOMATIC SHUTDOWN 10S ± 1S	BATTERY AAA X 2	ACCURACY(CLINICAL) 0.3°C / 0.6°F
MEASURING RANGE		
FOREHEAD MODE	22 - 43 °C / 71.6 - 109.4 °F	
EAR MODE	34 - 43 °C / 93.2 - 109.4 °F	
OBJECT MODE	0 - 100 °C / 32.0 - 212.0 °F	
Forehead/Ear Auto Switch Magnetic Probe Cap Child Mode 10s Auto Shut Off Smart Fever Alert Silent Mode		

NEW ARRIVAL



JPD-FR204

DISPLAY HD LED	MEMORY STORAGE 20 SETS DATA	TEMPERATURE UNITS °C/°F
AUTOMATIC SHUTDOWN 10S ± 1S	BATTERY AAA X 2	ACCURACY(CLINICAL) 0.3°C / 0.6°F
MEASURING RANGE		
FOREHEAD MODE	22 - 43 °C / 71.6 - 109.4 °F	
OBJECT MODE	0 - 100 °C / 32.0 - 212.0 °F	
 Child Mode 10s Auto Shut Off Smart Fever Alert Silent Mode		



JPD-FR410

DISPLAY HD LCD	MEMORY STORAGE 20 SETS DATA	TEMPERATURE UNITS °C/°F
AUTOMATIC SHUTDOWN 10S ± 1S	BATTERY AAA X 2	ACCURACY(CLINICAL) 0.3°C / 0.6°F
MEASURING RANGE		
FOREHEAD MODE	22 - 43 °C / 71.6 - 109.4 °F	
EAR MODE	34 - 43 °C / 93.2 - 109.4 °F	
OBJECT MODE	0 - 100 °C / 32.0 - 212.0 °F	
Child Mode 10s Auto Shut Off Smart Fever Alert Silent Mode		





JPD-FR409-BT

DISPLAY HD LCD	MEMORY STORAGE 20 SETS DATA	TEMPERATURE UNITS °C/°F
AUTOMATIC SHUTDOWN 10S ± 1S	BATTERY AAA X 2	ACCURACY(CLINICAL) 0.3°C / 0.6°F
MEASURING RANGE		
FOREHEAD MODE	22 - 43 °C / 71.6 - 109.4 °F	
EAR MODE	34 - 43 °C / 93.2 - 109.4 °F	
OBJECT MODE	0 - 100 °C / 32.0 - 212.0 °F	
Magnetic Probe Cap Child Mode 10s Auto Shut Off Smart Fever Alert Silent Mode		



JPD-FR409

DISPLAY HD LCD	MEMORY STORAGE 20 SETS DATA	TEMPERATURE UNITS °C/°F
AUTOMATIC SHUTDOWN 10S ± 1S	BATTERY AAA X 2	ACCURACY(CLINICAL) 0.3°C / 0.6°F
MEASURING RANGE		
FOREHEAD MODE	22 - 43 °C / 71.6 - 109.4 °F	
EAR MODE	34 - 43 °C / 93.2 - 109.4 °F	
OBJECT MODE	0 - 100 °C / 32.0 - 212.0 °F	
Magnetic Probe Cap Child Mode 10s Auto Shut Off Smart Fever Alert Silent Mode		





JPD-FR400

DISPLAY HD LCD	MEMORY STORAGE 20 SETS DATA	TEMPERATURE UNITS °C/°F
AUTOMATIC SHUTDOWN 10S ± 1S	BATTERY AAA X 2	ACCURACY(CLINICAL) 0.3°C / 0.6°F
MEASURING RANGE		
FOREHEAD/EAR MODE	32 - 42.2 °C / 89.6 - 107.9 °F	
OBJECT MODE	0.0 - 100 °C / 32.0 - 212.0 °F	
1S 1 Second Reading 10s Auto Shut Off Smart Fever Alert		



JPD-FR300

DISPLAY HD LCD	MEMORY STORAGE 20 SETS DATA	TEMPERATURE UNITS °C/°F
AUTOMATIC SHUTDOWN 10S ± 1S	BATTERY AAA X 2	ACCURACY(CLINICAL) 0.3°C / 0.6°F
MEASURING RANGE		
FOREHEAD MODE	22 - 43 °C / 71.6 - 109.4 °F	
EAR MODE	34 - 43 °C / 93.2 - 109.4 °F	
OBJECT MODE	0 - 100 °C / 32.0 - 212.0 °F	
 Forehead/Ear Auto Switch	 Magnetic Probe Cap	 Child Mode
 10s Auto Shut Off	 Smart Fever Alert	 Silent Mode



JPD-FR100+

DISPLAY HD LCD	MEMORY STORAGE 20 SETS DATA	TEMPERATURE UNITS °C/°F
AUTOMATIC SHUTDOWN 10S ± 1S	BATTERY AAA X 2	ACCURACY(CLINICAL) 0.3°C / 0.6°F
MEASURING RANGE		
FOREHEAD/EAR MODE	32 - 42.2 °C / 89.6 - 107.9 °F	
OBJECT MODE	0.0 - 100 °C / 32.0 - 212.0 °F	
 10s Auto Shut Off	 Smart Fever Alert	 1S Reading

NEW ARRIVAL



JPD-FR415

DISPLAY HD LCD	MEMORY STORAGE 20 SETS DATA	TEMPERATURE UNITS °C/°F
AUTOMATIC SHUTDOWN 10S ± 1S	BATTERY AAA X 2	ACCURACY(CLINICAL) 0.3°C / 0.6°F
MEASURING RANGE		
EAR MODE	34 - 43 °C / 93.2 - 109.4 °F	
 Magnetic Probe Cap	 Child Mode	 10s Auto Shut Off
 Smart Fever Alert	 Silent Mode	

NEW ARRIVAL



JPD-FR205

DISPLAY HD LCD	MEMORY STORAGE 20 SETS DATA	TEMPERATURE UNITS °C/°F
AUTOMATIC SHUTDOWN 10S ± 1S	BATTERY AAA X 2	ACCURACY(CLINICAL) 0.3°C / 0.6°F
MEASURING RANGE		
FOREHEAD MODE	22 - 43 °C / 71.6 - 109.4 °F	
OBJECT MODE	0 - 100 °C / 32.0 - 212.0 °F	
 Distance Detection	 Magnetic Battery Cap	 Hygienic Non-Touch
 10s Auto Shut Off	 Smart Fever Alert	 Silent Mode



JPD-FR401

DISPLAY HD LCD	MEMORY STORAGE 20 SETS DATA	TEMPERATURE UNITS °C/°F
AUTOMATIC SHUTDOWN 10S ± 1S	BATTERY AAA X 2	ACCURACY(CLINICAL) 0.3°C / 0.6°F
MEASURING RANGE		
FOREHEAD/EAR MODE	32 - 42.2 °C / 89.6 - 107.9 °F	
OBJECT MODE	0.0 - 100 °C / 32.0 - 212.0 °F	
 1S Reading	 10s Auto Shut Off	 Smart Fever Alert

NEW ARRIVAL



FTA10

DISPLAY HD LCD	MEMORY STORAGE 20 SETS DATA	TEMPERATURE UNITS °C/°F
AUTOMATIC SHUTDOWN 10S ± 1S	BATTERY AAA X 2	ACCURACY(CLINICAL) 0.3°C / 0.6°F
MEASURING RANGE		
FOREHEAD MODE	22 - 43 °C / 71.6 - 109.4 °F	
OBJECT MODE	0 - 100 °C / 32.0 - 212.0 °F	
 Hygienic Non-Touch	 Child Mode	 10s Auto Shut Off
 Smart Fever Alert	 Silent Mode	



JPD-FR203

DISPLAY HD LCD	MEMORY STORAGE 20 SETS DATA	TEMPERATURE UNITS °C/°F
AUTOMATIC SHUTDOWN 10S ± 1S	BATTERY AAA X 2	ACCURACY (CLINICAL) 0.3°C / 0.6°F
MEASURING RANGE		
FOREHEAD MODE	22 - 43 °C / 71.6 - 109.4 °F	
OBJECT MODE	0 - 100 °C / 32.0 - 212.0 °F	
Hygienic Non-Touch	Child Mode	10s Auto Shut Off
	Smart Fever Alert	Silent Mode



JPD-FR202

DISPLAY HD LCD	MEMORY STORAGE 20 SETS DATA	TEMPERATURE UNITS °C/°F
AUTOMATIC SHUTDOWN 10S ± 1S	BATTERY AAA X 2	ACCURACY (CLINICAL) 0.3°C / 0.6°F
MEASURING RANGE		
FOREHEAD MODE	22 - 43 °C / 71.6 - 109.4 °F	
OBJECT MODE	0 - 100 °C / 32.0 - 212.0 °F	
Hygienic Non-Touch	10s Auto Shut Off	Smart Fever Alert
	Silent Mode	

Infrared Thermometer Parameters											
Product Model	JPD-FR400	JPD-FR409-BT	JPD-FR403	JPD-FR302	JPD-FR301						
Display Screen	HD LCD				HD LED						
Forehead Mode	32.0-42.2 °C 89.6-107.9 °F	22.0-43.0 °C / 71.6-109.4 °F									
Ear Mode	0.0-100.0 °C 32.0-212.0 °F	34.0-43.0 °C / 93.2-109.4 °F									
Object Mode	0.0-100.0 °C / 32.0-212.0 °F										
Child Mode	—										
Silent Mode	—										
Backlight/ Indicator					—						
Clinical Accuracy	± 0.3 °C / ± 0.6 °F										
Auto Shut Off	10s ± 1s										
Features	❶ Fever warning ❷ °C/°F adjustable ❸ Memory storage: 20 sets data ❹ Batteries: AAA × 2 ❺ Forehead mode: for all ages ❻ Ear mode for children: ≥ 3 months old ❼ Child mode: ≤ 12 years old										

