



Official Website



JUMPER Health

Shenzhen Jumper Medical Equipment Co., Ltd.

Address: D Building, No. 71, Xintian Road, Fuyong Street, Baoan, Shenzhen, Guangdong, China

Tel: +86 755-26696279 **E-mail:** info@jumper-medical.com **Web:** <https://jumpermed.com>

PRODUCT BROCHURE

SINCE
1998



JUMPER

Healthcare within mind

<https://jumpermed.com>

TENS THERAPY DEVICE



- Bluetooth
- APP Control
- 5 Modes 15 Levels
- One Button Control
- Automatic Shutdown
- Relieve Pain
- Rechargeable Battery


RELIEVE PAIN, SOOTHE MIND.

Small, portable and powerful TENS devices from JUMPER offer drug-free pain relief for sore or aching muscles of the lower back, shoulder, arm, leg and more.

The unit also features one or two independent stimulation channels and four self-adhesive electrodes. It offers a wide range of functions for increasing general well-being, pain relief, maintaining physical fitness, relaxation, muscle revitalization and combating tiredness. For these purposes, you can either choose from pre-set programs or specify your own to suit your individual needs.

FEATURES

- Pain relief
- Pad off Alert
- Automatic shutdown
- Low battery indication
- Ultra-thin & lightweight
- Easy to charge and use
- Various massage methods
- Powerful & effective program







JPD-ES100

TREATMENT TIME	PULSE FREQUENCY	INTENSITY
1-60 mins	1-120HZ	0-15

MAX OUTPUT INTENSITY	POWER SUPPLY TYPE	PULSE WIDTH
30V at 500Ω	LITHIUM	0.15-0.20ms

DIMENSION(mm)	MAIN UNIT WEIGHT	
MAIN UNIT	44 x 44 x 13	50g
ELECTRODE PAD	φ90 / φ60	



APP ControlOne Button ControlPad Off Alert5 ModesRechargeable BatteryLow Battery Indication



JPD-ES200

TREATMENT TIME	PULSE FREQUENCY	INTENSITY
5-60 mins	1-120HZ	0-15

MAX OUTPUT INTENSITY	POWER SUPPLY TYPE	PULSE WIDTH
48mA at 500Ω	DRY BATTERIES	0.15-0.20ms

DIMENSION(mm)	MAIN UNIT WEIGHT	
MAIN UNIT	121.8 x 63.6 x 22.2	86g (WITHOUT BATTERY)
ELECTRODE PAD	50 x 50	



5 Modes15 Levels IntensityPad Off AlertAvailable For 2 UsersLow Battery Indication

TENSEMS

MESSAGE

EMS



JPD-ES220

TREATMENT TIME	PULSE FREQUENCY	INTENSITY
5-60 mins	1-100HZ ± 10%	0-20

MAX OUTPUT INTENSITY	POWER SUPPLY TYPE	PULSE WIDTH
50V at 500Ω	LITHIUM	80μs

DIMENSION(mm)	MAIN UNIT WEIGHT	
MAIN UNIT	138x 68.7 x 13	123.9g
ELECTRODE PAD	50 x 50	




Type-CRechargeable Battery20 Modes20 Levels IntensityAvailable For 2 UsersPad Off Alert



A


JPD-ES220

TREATMENT TIME	PULSE FREQUENCY	INTENSITY
5-60 mins	1-100HZ ± 10%	0-20
MAX OUTPUT INTENSITY	POWER SUPPLY TYPE	PULSE WIDTH
50V at 500Ω	LITHIUM	80μs
DIMENSION(mm)	MAIN UNIT WEIGHT	
MAIN UNIT	138x 68.7 x 13	
ELECTRODE PAD	50 x 50	



JPD-ES230

TREATMENT TIME	PULSE FREQUENCY	INTENSITY
5-60 mins	2-120HZ	0-40
MAX OUTPUT INTENSITY	POWER SUPPLY TYPE	PULSE WIDTH
50V at 500Ω	DRY BATTERIES	0.05-0.30ms
DIMENSION(mm)	MAIN UNIT WEIGHT	
MAIN UNIT	130 x 64 x 23	
ELECTRODE PAD	50 x 50	
98g (WITHOUT BATTERY)		



A

JPD-ES210

TREATMENT TIME	PULSE FREQUENCY	INTENSITY
5-60 mins	1-120HZ	0-15
MAX OUTPUT INTENSITY	POWER SUPPLY TYPE	PULSE WIDTH
48mA at 500Ω	DRY BATTERIES	0.1-0.2ms
DIMENSION(mm)	MAIN UNIT WEIGHT	
MAIN UNIT	132.3 x 64.8 x 19.6	
ELECTRODE PAD	98 x 64	
	105g (WITHOUT BATTERY)	



LPA11


TREATMENT TIME	PULSE FREQUENCY	INTENSITY
5-60 mins	2-100HZ	0-40
MAX OUTPUT INTENSITY	POWER SUPPLY TYPE	PULSE WIDTH
50V at 500Ω	DRY BATTERIES	0.05-0.30ms
DIMENSION(mm)	MAIN UNIT WEIGHT	
MAIN UNIT	113x 62 x 22.5	
ELECTRODE PAD	50 x 50	



B


JPD-ES210

TREATMENT TIME	PULSE FREQUENCY	INTENSITY
5-60 mins	2-120HZ	0-15
MAX OUTPUT INTENSITY	POWER SUPPLY TYPE	PULSE WIDTH
48mA at 500Ω	DRY BATTERIES	0.05-0.30ms
DIMENSION(mm)	MAIN UNIT WEIGHT	
MAIN UNIT	132.3 x 64.8 x 19.6	
ELECTRODE PAD	98 x 64	
	105g (WITHOUT BATTERY)	



LPA12

TREATMENT TIME	PULSE FREQUENCY	INTENSITY
5-60 mins	2-120HZ	0-50
MAX OUTPUT INTENSITY	POWER SUPPLY TYPE	PULSE WIDTH
50V at 500Ω	LITHIUM	0.05-0.30ms
DIMENSION(mm)	MAIN UNIT WEIGHT	
MAIN UNIT	151 x 78 x 32	
ELECTRODE PAD	50 x 50	
	147g (WITHOUT BATTERY)	



★NEW★
ARRIVAL

TENS

MASSAGE

EMS

FITNESS



LPB12

TREATMENT TIME	PULSE FREQUENCY	INTENSITY
5-60 mins	2-120HZ	0-50
MAX OUTPUT INTENSITY	POWER SUPPLY TYPE	PULSE WIDTH
50V at 500Ω	LITHIUM	0.05-0.30ms
DIMENSION(mm)	MAIN UNIT WEIGHT	
MAIN UNIT	151 x 82 x 33	
ELECTRODE PAD	50 x 50	
128g (WITHOUT BATTERY)		

79 Programs

50 Levels Intensity

Type-C

Pad Off Alert

Available For 2 Users

Low Battery Indication

TENS Therapy Device Parameters	 B  A  A  				
 Product Model	JPD-ES220	JPD-ES220	JPD-ES210	JPD-ES200	JPD-ES100
 TENS EMS MESSAGE	TENS - EMS - MESSAGE	TENS			
 Treatment Time	5-60 mins		5-60 mins	5-60 mins	1-60 mins
 Power Supply Type	Lithium		Dry batteries	Dry batteries	Lithium
 Max Output Intensity (at 500Ω)	50 V		48 mA	48 mA	30 V
 Pulse Width	80μs		0.1-0.2ms	0.15-0.20ms	
 Pulse Frequency	1-100Hz		1-120Hz		
 Low Battery Indication	✓				
 Number of Channels	2				1
 Auto Power Off	✓				
 Bluetooth	—				
 Main Unit Size (mm)	138 x 68.7 x 13		132.3 x 64.8 x 19.6	121.8 x 63.6 x 22.2	44 x 44 x 13
 Electrode Pad Size (mm)	50 x 50		98 x 64	50 x 50	φ90 / φ60
 Main Unit Weight (g)	123.9		105	86	50
 Modes	25	20	14	5	

TENS Therapy Device Parameters		B				
Product Model	JPD-ES210	JPD-ES230	LPA11	LPA12	LPB12	
TENS EMS MASSAGE	TENS - EMS - MASSAGE			TENS - EMS - MASSAGE - FITNESS		
Treatment Time	5-60 mins					
Power Supply Type	Dry batteries			Lithium		
Max Output Intensity (at 500Ω)	48 mA	50 V				
Pulse Width	0.05-0.3ms					
Pulse Frequency	2-120Hz					
Low Battery Indication	✓					
Number of Channels	2					
Auto Power Off	✓					
Bluetooth	—					
Main Unit Size (mm)	132.3 x 64.8 x 19.6	130 x 64 x 23	113 x 61 x 22	151 x 78 x 32	151 x 82 x 33	
Electrode Pad Size (mm)	98 x 64	50 x 50	50 x 50	50 x 50	50 x 50	
Main Unit Weight (g)	105	98	95	147	128	
Modes	41	25	25	79	79	