



Official Website



JUMPER Health

Shenzhen Jumper Medical Equipment Co., Ltd.

Address: D Building, No. 71, Xintian Road, Fuyong Street, Baoan, Shenzhen, Guangdong, China

Tel: +86 755-26696279 **E-mail:** info@jumper-medical.com **Web:** <https://jumpermed.com>

PRODUCT BROCHURE

SINCE
1998



JUMPER

Healthcare within mind

<https://jumpermed.com>

MESH NEBULIZER

COMPRESSOR

NEBULIZER



- 
Lightweight
and Portable
- 
Well
absorbed
- 
High
efficiency
- 
Low noise
- 
Independent
nebulizing
components

COMFORTABLE NEBULIZATION

Our nebulizer is high-performance and convenient. It fits various types of nebulized medication and can nebulize large particles of medication and viscous liquids. The nebulizer offers an affordable and efficient alternative for patients experiencing asthma, COPD or other respiratory disease conditions.



CTB10

NEW
ARRIVAL



MTA10

POWER	MEDICATION CAPACITY	NEBULIZATION RATE
DC 3V	14 ML MAXIMUM	≥ 0.2 ML/MIN

SOUND LEVEL	PARTICLE SIZE
≤ 45 DB	1-5µm≥ 50 % (MMAD:3.5µm)

DIMENSION(mm)	Net Weight
MAIN UNIT	104g

- 
Lightweight
and Portable
- 
One-button
Operation
- 
Type-C / Battery
Power Supply
- 
Low noise
- 
14 ml Large
Capacity
- 
Self-Cleaning

NEW
ARRIVAL



CTA10

POWER	MEDICATION CAPACITY	NEBULIZATION RATE
AC 220V	8 ML MAXIMUM	≥ 0.2 ML/MIN

SOUND LEVEL	PARTICLE SIZE
≤ 50 DB	3.5 µm ±25%

DIMENSION(mm)	Net Weight
MAIN UNIT	1450g

- 
Easy to use
- 
Well
absorbed
- 
Adjustable
atomization
- 
Low noise
- 
Effective
and fast
medication
delivery
- 
More fully

NEW
ARRIVAL



CTB10

POWER	MEDICATION CAPACITY	NEBULIZATION RATE
AC 220V	8 ML MAXIMUM	≥ 0.2 ML/MIN

SOUND LEVEL	PARTICLE SIZE
≤ 50 DB	3.0 µm ±10%

DIMENSION(mm)	Net Weight
MAIN UNIT	170g

- 
Lightweight
and portable
- 
Well
absorbed
- 
Adjustable
atomization
- 
Low noise
- 
Effective
and fast
medication
delivery
- 
More fully

Nebulizer Comparison Table Parameters		  		
Product Model	MTA10	CTA10	CTB10	
Product Type	Mesh Nebulizer	Compressor Nebulizer	Compressor Nebulizer	
Nebulization Rate	≥ 0.2ml/min	≥ 0.2ml/min	≥ 0.2ml/min	
Power Supply	DC 5V / AAx2	AC 220V/50Hz	AC 100-240V-50 / 60Hz 0.4A	
Particles	1-5μm ≥ 50% (MMAD: 3.5μm)	3.5μm±25%	3.0μm±10%	
Sound	≤ 45dB	≤ 50 dB	≤ 50 dB	
Dimension (mm)	62.1 x 45 x 125.5	167 x 167 x 102	100 x 100 x 49.5	
Net Weight (g)	104	1450	170	
Accessories	Nebulization cup Adult mask Child mask Mouthpiece Headstrap	Nebulization cup Adult mask Child mask Mouthpiece Headstrap Airtube Air filters	Nebulization cup Adult mask Child mask Mouthpiece Headstrap Airtube Air filters	
Applicable Users	Elderly Adults Children	Elderly Adults Children	Elderly Adults Children	
Applicable Scenarios	Family	Hospital Clinic Family	Hospital Clinic Family	