



Official Website



JUMPER Health

Shenzhen Jumper Medical Equipment Co., Ltd.

Address: D Building, No. 71, Xintian Road, Fuyong Street, Baoan, Shenzhen, Guangdong, China

Tel: +86 755-26696279 **E-mail:** info@jumper-medical.com **Web:** <https://jumpermed.com>

PRODUCT BROCHURE

SINCE
1998



JUMPER

Healthcare within mind

<https://jumpermed.com>

WEIGHT SCALE

BODY FAT SCALE

- 

BMI
- 

Body Fat Rate
- 

Body Fat
- 

Body Water
- 

Automatic Power-Off

PROFESSIONAL MEASUREMENT, SCIENTIFICALLY MANAGE HEALTH.

These products are designed to measure body weight and body fat rate.

It can be connected to a smart phone via Bluetooth, and the measured data will be transmitted to the phone synchronously, which means that users can check their weight anytime by turning on the APP.

These products help users to develop good habits of diet management and exercise for the purpose of body management.

FEATURES


- Extra-large platform.

•Elegant and stylish design.

•Extra-high weight capacity.
- Accurate weight measurement.

•Durable tempered glass surface.

•Large easy-to-read LED/LCD display.
- The data of measurement can be checked on APP (JUMPER-HEALTH) via Bluetooth.







JPD-BFS710

DISPLAY	DIMENSIONS	POWER SUPPLY	
HD LED	325 x 325 x 23.5 mm	AAA x 4	

NET WEIGHT	MAX CAPACITY	DIVISION	UNIT
3 KG	180 KG / 396 LB	0.1 KG / 0.2 LB	KG / LB / ST

FUNCTION

LOW POWER & OVERLOAD INDICATION

AUTO POWER OFF IN 10 SECONDS



BMI



Body Fat Rate



Body Fat



Body Water



Automatic Power-Off





JPD-BFS100

DISPLAY	DIMENSIONS	POWER SUPPLY	UNIT
HD LED	372 x 298 x 29 mm	AAA x 3	KG / LB

NET WEIGHT	MAX CAPACITY	DIVISION
3 KG	180 KG / 396 LB	0.1 KG / 0.2 LB

FUNCTION

LOW POWER & OVERLOAD INDICATION

AUTO POWER OFF IN 10 SECONDS

BMI



Body Fat Rate



Body Fat



Body Water



Automatic Power-Off







JPD-BFS200B



JPD-BFS200D

JPD-BFS200A

DISPLAY	DIMENSIONS	POWER SUPPLY	
HD LED	280 x 280 x 23.5 mm	AAA x 3	

NET WEIGHT	MAX CAPACITY	DIVISION	UNIT
1.3 KG	200 KG / 440 LB	0.1 KG / 0.2 LB	KG / LB

FUNCTION

LOW POWER & OVERLOAD INDICATION

AUTO POWER OFF IN 10 SECONDS

BMI



Body Fat Rate



Body Fat



Body Water



Automatic Power-Off



JPD-BFS300



DISPLAY	DIMENSIONS	UNIT	POWER SUPPLY
HD LCD	280 x 280 x 23.5 mm	KG / LB / ST	AAA x 4

NET WEIGHT	MAX CAPACITY	DIVISION
1.3 KG	200 KG / 440 LB	0.1 KG / 0.2 LB

FUNCTION
LOW POWER & OVERLOAD INDICATION
AUTO POWER OFF IN 10 SECONDS

BMI

Body Fat Rate

Body Fat

Body Water

Automatic Power-Off



JPD-700A



DISPLAY	DIMENSIONS	MAX CAPACITY
HD LED	375 x 300 x 24.5 mm	200 KG / 440 LB

NET WEIGHT	DIVISION	UNIT	POWER SUPPLY
2.1 KG	0.1 KG / 0.2 LB	KG / LB	AAA x 3

FUNCTION
LOW POWER & OVERLOAD INDICATION
AUTO POWER OFF IN 10 SECONDS

Low Battery Indication

Unit Conversion

Automatic Power-Off

LED Screen Display

Top Quality Material

A



JPD-BS300



DISPLAY	DIMENSIONS	MAX CAPACITY
HD LCD	280 x 280 x 23.5 mm	200 KG / 440 LB / 31 ST

NET WEIGHT	DIVISION	UNIT	POWER SUPPLY
1.3 KG	0.1 KG / 0.2 LB	KG / LB / ST	AAA x 3

FUNCTION
LOW POWER & OVERLOAD INDICATION
AUTO POWER OFF IN 10 SECONDS

Low Battery Indication

Unit Conversion

Automatic Power-Off

LED Screen Display

Top Quality Material



JPD-700A



DISPLAY	DIMENSIONS	MAX CAPACITY
HD LED	375 x 300 x 24.5 mm	200 KG / 440 LB

NET WEIGHT	DIVISION	UNIT
2.1 KG	0.1 KG / 0.2 LB	KG / LB

FUNCTION
LOW POWER & OVERLOAD INDICATION
AUTO POWER OFF IN 10 SECONDS

Low Battery Indication

Unit Conversion

Automatic Power-Off

LED Screen Display

Top Quality Material

B



JPD-BS200



JPD-BS201

DISPLAY	DIMENSIONS	MAX CAPACITY
HD LED	280 x 280 x 23.5 mm	200 KG / 440 LB

NET WEIGHT	DIVISION	UNIT	POWER SUPPLY
1.3 KG	0.1 KG / 0.2 LB	KG / LB	AAA x 3

FUNCTION
LOW POWER & OVERLOAD INDICATION
AUTO POWER OFF IN 10 SECONDS

LED Screen Display

Top Quality Material

Automatic Power-Off

JPD-700C



DISPLAY	DIMENSIONS	MAX CAPACITY
HD LED	315 x 283 x 24.5 mm	180 KG / 396 LB

NET WEIGHT	DIVISION	UNIT	POWER SUPPLY
1.6 KG	0.1 KG / 0.2 LB	KG / LB	AAA x 4

FUNCTION
LOW POWER & OVERLOAD INDICATION
AUTO POWER OFF IN 10 SECONDS

Low Battery Indication

Unit Conversion

LED Screen Display

Top Quality Material

Automatic Power-Off

Weight & Body Fat Scale Parameters												
	Product Model	JPD-BFS710	JPD-BFS100	JPD-BFS200A/B/D	JPD-BFS300	JPD-BS300	JPD-BS200/201		JPD-700A	JPD-700A	JPD-700C	
	Display	HD LED		HD LED	HD LCD	HD LCD	HD LED		HD LED			
	Bluetooth	✖							✖		—	
	Power Supply	AAA x 4	AAA x 3	AAA x 3	AAA x 4	AAA x 3			AAA x 3		AAA x 4	
Dimensions (mm)	325 x 325 x 23.5	372 x 298 x 29	280 x 280 x 23.5					375x300x24.5		315 x 283 x 24.5		
Net Weight (Kg)	3	3	1.3					2.4		1.6		
Auto Power Off	10s							10s				
Max Capacity	180 KG / 396 LB		200 KG / 440 LB					200 KG / 440 LB		180 KG / 396 LB		
Indication	LOW POWER & OVERLOAD							LOW POWER & OVERLOAD				
Back Side												