



Official Website



JUMPER Health

Shenzhen Jumper Medical Equipment Co., Ltd.

Address: D Building, No. 71, Xintian Road, Fuyong Street, Baoan, Shenzhen, Guangdong, China

Tel: +86 755-26696279 **E-mail:** info@jumper-medical.com **Web:** <https://jumpermed.com>

PRODUCT BROCHURE

SINCE
1998



JUMPER

Healthcare within mind

<https://jumpermed.com>

About Us

Founded in **1998**, Shenzhen Jumper Medical Equipment Co., Ltd, is a national high-tech enterprise, integrating research and development, production, sales, and service of intelligent medical equipment. As a professional manufacturer, we provide homecare devices & medical equipment and solutions for maternal and child healthcare.

MISSION

Technology Drives Maternal and Child Health.

PROFESSIONAL STRENGTH

Completed R&D, production, marketing and service network.

- **3** R&D centers: Shenzhen, Beijing, Changsha.
- **14** subsidiaries (branches) across the country.
- Over **16,000** square meters of office area.

MARKET LAYOUT

- More than **160** countries and regions.
- AngelSounds series: **30%** European and US markets.

CORE VALUE

Integrity, Dedication, Simplicity, Efficiency, Customer first, Embrace change.

MEDICAL EQUIPMENT

- Fetal & Maternal Monitor
- Multi-parameter Patient Monitor
- Fetal Monitor
- AngelSounds Fetal Doppler
- Pocket Fetal Doppler
- Infrared Thermometer
- Fingertip Pulse Oximeter
- Blood Pressure Monitor
- Weight Scale & Body Fat Scale
- Tens Therapy Device
- Mesh Nebulizer & Compressor Nebulizer
- Pelvic Floor Rehabilitation Device



International Certificates

PRODUCT

EUROPE



US



CANADA



CN



BRAZIL



WIRELESS

EUROPE



US



CANADA



ISO




PACKAGE



BUSINESS SOCIAL COMPLIANCE INITIATIVE



Development




Contents

	ANGELSOUNDS FETAL DOPPLER Page 01-04	01		POCKET FETAL DOPPLER Page 05-13	02
	FINGERTIP PULSE OXIMETER Page 14-23	03		INFRARED THERMOMETER Page 24-33	04
	BLOOD PRESSURE MONITOR Page 34-44	05		TENS THERAPY DEVICE Page 45-51	06
	PELVIC FLOOR REHABILITATION DEVICE Page 52-54	07		MESH NEBULIZER & COMPRESSOR NEBULIZER Page 55-57	08
	WEIGHT SCALE & BODY FAT SCALE Page 58-63	09		REMOTE PATIENT MONITORING DEVICES Page 64-68	10

AngelSounds®

FETAL DOPPLER



CATCH BABY’S HEARTBEAT.

AngelSounds Fetal Doppler is a good helper for mothers to perceive the baby. With the function of monitoring the fetal heart rate, it's also a baby's health guardian. Sixteen weeks into your pregnancy, the sounds of heartbeats, hiccups and kicks can now be heard. Whether you face a high risk of pregnancy or just want to ensure that your angel is safe and well, it's helpful to use Angelsounds Fetal Doppler.




JPD-100S9

ULTRASOUND FREQUENCY	BATTERY	DATA TRANSMISSION
2.5 MHz	IEC6F22 9V Alkaline	Wired (Alternative) Bluetooth

WEIGHT	DIMENSIONS(mm)	MEASURING MODE
155g	119 x 65 x 54	REAL-TIME

FUNCTION
ADJUSTABLE VOLUME




JPD-100S+

ULTRASOUND FREQUENCY	BATTERY	DATA TRANSMISSION
2.0 MHz	3.7V Lithium (RECHARGEABLE)	Bluetooth

NET WEIGHT	DIMENSIONS(mm)	MEASURING MODE
99.5g	80 x 80 x 50	REAL-TIME

FUNCTION
ADJUSTABLE VOLUME
APP CONNECTION




JPD-100S(Mini)

ULTRASOUND FREQUENCY	BATTERY	DATA TRANSMISSION
3.0 MHz	IEC6F22 9V Alkaline	Wired

NET WEIGHT	DIMENSIONS(mm)	MEASURING MODE
69.5g (WITHOUT BATTERY)	104.5 x 50 x 70	REAL-TIME

FUNCTION
ADJUSTABLE VOLUME





JPD-100S

ULTRASOUND FREQUENCY 3.0 MHz	BATTERY IEC6F22 9V Alkaline	DATA TRANSMISSION Wired
WEIGHT 140g	DIMENSIONS(mm) 104 x 120 x 61	MEASURING MODE REAL-TIME
FUNCTION ADJUSTABLE VOLUME		


FHR


Clinically Tested


From 16th Week


Battery Included





JPD-200S

ULTRASOUND FREQUENCY 3.0 MHz	BATTERY IEC6F22 9V Alkaline	DATA TRANSMISSION Wired
WEIGHT 150g	DIMENSIONS(mm) 100 x 85 x 50	MEASURING MODE REAL-TIME
FUNCTION ADJUSTABLE VOLUME		


FHR


Clinically Tested


From 16th Week


Battery Included


Wired Connection

AngelSounds Fetal Doppler Parameters					
 Product Model	JPD-100S+	JPD-100S	JPD-100S (mini)	JPD-200S	JPD-100S9
 Ultrasound Frequency	2.0 MHz	3.0 MHz			2.5 MHz
 Display Screen	—				
 FHR Value	—				
 Volume Display	—				
 Volume Adjustable	✓				
 Battery Status	—				
 Data Transmission		—	Wired		Wired / Bluetooth
 Measuring Mode	 Real-time				
 Backlight	—				
 Auto Power Off	✓	—			
 Battery Type	3.7V Lithium (RECHARGEABLE)	IEC6F22 9V ALKALINE			
 Dimensions (mm)	80 x 80 x 50	104 x 120 x 61	104.5 x 50 x 70	100 x 85 x 50	119 x 65 x 54
 Weight(g) (Include Batteries)	99.5	140	69.5 (WITHOUT BATTERY)	150	155
 Features	 FHR Measurement Range: 50-210 ± 2 bpm				

POCKET FETAL DOPPLER



From
9th Week



FHR




Battery
Included



TFT
Screen



1 min Auto
Power Off



HEAR THE BABY'S HEARTBEAT.

Pocket Fetal Doppler is a hand-held fetal heart rate detecting device.

It displays fetal heart rate on a digital LCD screen /TFT screen with easy and convenient operation, which can be used in hospitals, clinics and home for daily self-check.



SHB10

DISPLAY LED	ULTRASOUND FREQUENCY 3.0 MHz	BATTERY AAA X 5
-----------------------	--	---------------------------

NET WEIGHT 156.5g (WITHOUT BATTERY)	DIMENSIONS(mm) 154 x 69 x 32	FHR MEASUREMENT RANGE 50-240 BPM
--	--	--

PERFORMANCE BATTERY STATUS FHR VALUE VOLUME	MEASURING MODE REAL-TIME	FUNCTION ADJUSTABLE VOLUME
---	------------------------------------	--------------------------------------



FHR



Clinically
Tested



From
16th Week



Battery
Included



LED
Screen



Adjustable
Volume



SHB20

DISPLAY LED	ULTRASOUND FREQUENCY 3.0 MHz	BATTERY AAA X 5
-----------------------	--	---------------------------

NET WEIGHT 156.5g (WITHOUT BATTERY)	DIMENSIONS(mm) 154 x 69 x 32	FHR MEASUREMENT RANGE 50-240 BPM
--	--	--

PERFORMANCE BATTERY STATUS FHR VALUE VOLUME	MEASURING MODE REAL-TIME	FUNCTION ADJUSTABLE VOLUME
---	------------------------------------	--------------------------------------



FHR



Clinically
Tested



From
16th Week




Battery
Included



LED
Screen



Adjustable
Volume



Sound
Alert Function



JPD-100F

DISPLAY TFT	ULTRASOUND FREQUENCY 2.5 MHz	BATTERY 3.7 V Lithium (RECHARGEABLE)
-----------------------	--	---

NET WEIGHT 250g	DIMENSIONS(mm) 139 x 82 x 28	FHR MEASUREMENT RANGE 50-210 ± 2 BPM
---------------------------	--	--

PERFORMANCE BATTERY STATUS WAVE FORM FHR VALUE VOLUME	MEASURING MODE REAL-TIME	FUNCTION ADJUSTABLE VOLUME
---	------------------------------------	--------------------------------------



FHR



Clinically
Tested



From
16th Week



Battery
Included



TFT
Screen



1 min Auto
Power Off



SHA10

DISPLAY TFT	ULTRASOUND FREQUENCY 3.0 MHz	BATTERY AA x 2
NET WEIGHT 130g	DIMENSIONS(mm) 142.6 x 40.5 x 42.5	FHR MEASUREMENT RANGE 50-210 ± 2 BPM
PERFORMANCE BATTERY STATUS WAVE FORM FHR VALUE VOLUME	MEASURING MODE REAL-TIME	FUNCTION ADJUSTABLE VOLUME BACKLIGHT
		
FHR	Clinically Tested	From 9th Week
		
Battery Included	TFT Screen	1 min Auto Power Off



JPD-100E2
(AAA X 3)

JPD-100E

DISPLAY TFT	ULTRASOUND FREQUENCY 2.5 MHz	BATTERY 3.7 V Lithium (RECHARGEABLE)
NET WEIGHT 250g (WITHOUT BATTERY)	DIMENSIONS(mm) 147 x 101.9 x 31.6	FHR MEASUREMENT RANGE 50-210 ± 2 BPM
PERFORMANCE BATTERY STATUS WAVE FORM FHR VALUE VOLUME	MEASURING MODE REAL-TIME	FUNCTION ADJUSTABLE VOLUME
		
FHR	Clinically Tested	From 16th Week
		
Battery Included	TFT Screen	1 min Auto Power Off



SHA20

DISPLAY TFT	ULTRASOUND FREQUENCY 3.0 MHz	BATTERY AA x 2
NET WEIGHT 130g	DIMENSIONS(mm) 142.6 x 40.5 x 42.5	FHR MEASUREMENT RANGE 50-210 ± 2 BPM
PERFORMANCE BATTERY STATUS WAVE FORM FHR VALUE VOLUME	MEASURING MODE REAL-TIME	FUNCTION ADJUSTABLE VOLUME BACKLIGHT
		
FHR	Clinically Tested	From 9th Week
		
Battery Included	TFT Screen	1 min Auto Power Off



JPD-100A2
(AA X 3)

JPD-100A

DISPLAY HD LCD	ULTRASOUND FREQUENCY 3.0 MHz	BATTERY 7.2V NI-MH (RECHARGEABLE)
NET WEIGHT 240g (WITHOUT BATTERY)	DIMENSIONS(mm) 135 x 83 x 33	FHR MEASUREMENT RANGE 50-210 ± 2 BPM
PERFORMANCE BATTERY STATUS FHR VALUE VOLUME	MEASURING MODE MANUAL COUNTING REAL-TIME AVERAGE	FUNCTION ADJUSTABLE VOLUME BACKLIGHT
		
FHR	Clinically Tested	From 16th Week
		
Battery Included	LCD Screen	1 min Auto Power Off



JPD-100E1
(AAA X 3)

JPD-100E+

DISPLAY TFT	ULTRASOUND FREQUENCY 2.5 MHz	BATTERY 3.7 V Lithium (RECHARGEABLE)
NET WEIGHT 250g	DIMENSIONS(mm) 147 x 101.9 x 31.6	FHR MEASUREMENT RANGE 50-210 ± 2 BPM
PERFORMANCE BATTERY STATUS WAVE FORM FHR VALUE VOLUME	MEASURING MODE REAL-TIME	FUNCTION ADJUSTABLE VOLUME
		
FHR	Clinically Tested	From 16th Week
		
Battery Included	TFT Screen	1 min Auto Power Off



JPD-100S6+

DISPLAY HD LCD	ULTRASOUND FREQUENCY 3.0 MHz	BATTERY AA x 2
NET WEIGHT 250g	DIMENSIONS(mm) 130 x 100 x 36	FHR MEASUREMENT RANGE 50-210 ± 2 BPM
PERFORMANCE BATTERY STATUS FHR VALUE	MEASURING MODE MANUAL COUNTING REAL-TIME AVERAGE	FUNCTION ADJUSTABLE VOLUME BACKLIGHT
		
FHR	Clinically Tested	From 16th Week
		
Battery Included	LCD Screen	1 min Auto Power Off



JPD-100S6-2

DISPLAY	ULTRASOUND FREQUENCY	BATTERY
TFT	3.0 MHz	AA x 2

NET WEIGHT	DIMENSIONS(mm)	FHR MEASUREMENT RANGE
250g	130 x 100 x 36	50-210 ± 2 BPM

PERFORMANCE	MEASURING MODE	FUNCTION
BATTERY STATUS WAVE FORM FHR VALUE	MANUAL COUNTING REAL-TIME AVERAGE	ADJUSTABLE VOLUME BACKLIGHT



FHR



Clinically
Tested



From
16th Week



Battery
Included



TFT
Screen



1 min Auto
Power Off



JPD-100B3
(AA X 2)

JPD-100B4

DISPLAY	ULTRASOUND FREQUENCY	BATTERY
TFT	2.5 MHz	9.6V NI-MH (RECHARGEABLE)

NET WEIGHT	DIMENSIONS(mm)	FHR MEASUREMENT RANGE
240g (WITHOUT BATTERY)	129 x 83 x 34	50-210 ± 2 BPM

PERFORMANCE	MEASURING MODE	FUNCTION
BATTERY STATUS WAVE FORM FHR VALUE VOLUME	REAL-TIME	ADJUSTABLE VOLUME BACKLIGHT



FHR



Clinically
Tested



From
16th Week



Battery
Included



TFT
Screen



1 min Auto
Power Off



JPD-100S6

DISPLAY	ULTRASOUND FREQUENCY	BATTERY
HD LCD	3.0 MHz	AA x 2

NET WEIGHT	DIMENSIONS(mm)	FHR MEASUREMENT RANGE
250g	130 x 100 x 36	50-210 ± 2 BPM

PERFORMANCE	MEASURING MODE	FUNCTION
BATTERY STATUS FHR VALUE	MANUAL COUNTING REAL-TIME AVERAGE	ADJUSTABLE VOLUME BACKLIGHT



FHR



Clinically
Tested



From
16th Week



Battery
Included



LCD
Screen



1 min Auto
Power Off



JPD-100B+

DISPLAY	ULTRASOUND FREQUENCY	BATTERY
HD LCD	2.5 MHz	7.4V Lithium (RECHARGEABLE)

NET WEIGHT	DIMENSIONS(mm)	FHR MEASUREMENT RANGE
250g (WITHOUT BATTERY)	195 x 190 x 52	50-210 ± 2 BPM

PERFORMANCE	MEASURING MODE	FUNCTION
BATTERY STATUS FHR VALUE	MANUAL COUNTING REAL-TIME AVERAGE	ADJUSTABLE VOLUME BACKLIGHT



FHR



Clinically
Tested



From
16th Week



Battery
Included



LCD
Screen



1 min Auto
Power Off



JPD-100B2
(AA X 2)

JPD-100B

DISPLAY	ULTRASOUND FREQUENCY	BATTERY
HD LCD	2.5 MHz	9.6V NI-MH (RECHARGEABLE)

NET WEIGHT	DIMENSIONS(mm)	FHR MEASUREMENT RANGE
240g (WITHOUT BATTERY)	129 x 83 x 34	50-210 ± 2 BPM

PERFORMANCE	MEASURING MODE	FUNCTION
BATTERY STATUS FHR VALUE	REAL-TIME	ADJUSTABLE VOLUME BACKLIGHT



FHR



Clinically
Tested



From
16th Week



Battery
Included



LCD
Screen



Adjustable
Volume

Pocket Fetal Doppler			
Parameters			
 Product Model	JPD-100B/JPD-100B2	JPD-100B3/JPD-100B4	
 Ultrasound Frequency	2.5 MHz		
 Display Screen	LCD	TFT	
 Bluetooth	—		
 Volume Display	—	✓	
 Wave Form	—	✓	
 Measuring Mode	① Real-time		
 Battery Type	9.6V NI-MH / AA x 2 (RECHARGEABLE)		
 Dimensions (mm)	129 x 83 x 34		
 Weight(g) (Include Batteries)	240 (WITHOUT BATTERY)		
 Features	① FHR Measurement Range: 50-210 ± 2 bpm ② Volume Adjustable ③ Battery Status ④ 1 min auto power off		

FINGERTIP PULSE OXIMETER



YOUR HEALTH WITHIN REACH

Fingertip pulse oximeter can measure the amount of oxygenated hemoglobin in your blood over a period of time.

It adopts Photoelectric Oxyhemoglobin Inspection Technology in accordance with Capacity Pulse Scanning Recording Technology, which can be used to measure human oxygen saturation and pulse rate through finger.

The device is suitable for daily care in family, hospitals, oxygen bars, community healthcare and physical care in sports, etc. It can be used before or after exercise, it is recommended to measure in a calm state.



FEATURES

- Multi-directional display reading.
- HD large screen provide clear measurement.
- High accuracy real-time measurement indication.
- Unique soft silicon housing provide comfortable measurement.
- With ALERT function will remind patients of parameter changes.
- Use bluetooth to transmit data to Android, IOS Devices & Windows Operation Systems.



SPA11

DISPLAY	PARAMETERS	APPLICATION RANGE	
TFT	PR / SpO ₂ / PI	ADULT / ADOLESCENT	

DIMENSION	DISPLAY DIRECTION	POWER SUPPLY	WEIGHT
62 x 33 x 32 mm	4 DIRECTIONS	AAA x 2	49.5g

RANGE & ACCURACY			INDICATION
SpO ₂	70 - 100%	± 2 %	REAL-TIME BATTERY STATUS VALUE LIMITED REMIND
PULSE RATE	25-250 BPM	± 2 BPM	




SPA10

DISPLAY	PARAMETERS	APPLICATION RANGE	
LED	PR / SpO ₂	ADULT / ADOLESCENT	

DIMENSION	POWER SUPPLY	WEIGHT
62 x 33 x 32 mm	AAA x 2	49.5g

RANGE & ACCURACY			INDICATION
SpO ₂	70 - 100%	± 2 %	LOW BATTERY REMIND VALUE LIMITED REMIND
PULSE RATE	25 - 250 BPM	± 2 BPM	




JPD-500I

DISPLAY	PARAMETERS	APPLICATION RANGE	
TFT	PR / SpO ₂ / PI	ADULT / ADOLESCENT	

DIMENSION	DISPLAY DIRECTION	POWER SUPPLY	WEIGHT
58 x 33.9 x 32 mm	4 DIRECTIONS	AAA x 2	49.5g

RANGE & ACCURACY			INDICATION
SpO ₂	70 - 100%	± 2 %	REAL-TIME BATTERY STATUS VALUE LIMITED REMIND
PULSE RATE	25-250 BPM	± 2 BPM	





Extra-Large
Screen



SPA20

DISPLAY	PARAMETERS	APPLICATION RANGE
LED	PR / SpO ₂	ADULT / ADOLESCENT

DIMENSION	DISPLAY DIRECTION	POWER SUPPLY	WEIGHT
67 x 35 x 31.7 mm	2 DIRECTIONS	AAA x 2	60g

RANGE & ACCURACY			INDICATION
SpO ₂	70 - 99 %	± 2 %	LOW BATTERY REMIND
PULSE RATE	25-250 BPM	± 2 BPM	VALUE LIMITED REMIND

Unique Finger Chamber Design

2 Display Directions

Pulse Rate Measuring

SpO₂ Measuring

Alert Function



Extra-Large
Screen



SPA30

DISPLAY	PARAMETERS	APPLICATION RANGE
LED	PR / SpO ₂	ADULT / ADOLESCENT

DIMENSION	DISPLAY DIRECTION	POWER SUPPLY	WEIGHT
62 x 33 x 31.5 mm	2 DIRECTIONS	AAA x 2	51.5g

RANGE & ACCURACY			INDICATION
SpO ₂	70 - 99 %	± 2 %	LOW BATTERY REMIND
PULSE RATE	25-250 BPM	± 2 BPM	VALUE LIMITED REMIND

Unique Finger Chamber Design

2 Display Directions

Pulse Rate Measuring

SpO₂ Measuring

Alert Function



JPD-500F

DISPLAY	PARAMETERS	APPLICATION RANGE
OLED	PR / SpO ₂ / PI	ADULT / ADOLESCENT

DIMENSION	DISPLAY DIRECTION	POWER SUPPLY	WEIGHT
62 x 37 x 32 mm	4 DIRECTIONS	AAA x 2	42.5g

RANGE & ACCURACY			INDICATION
SpO ₂	70 - 100%	± 2 %	REAL-TIME BATTERY STATUS
PULSE RATE	25-250 BPM	± 2 BPM	VALUE LIMITED REMIND

Multiple Display Modes

Pulse Rate Measuring

SpO₂ Measuring

Brightness Adjustable

Alert Function



Rotating Display
(Optional)



JPD-501F

DISPLAY	PARAMETERS	APPLICATION RANGE
OLED	PR / SpO ₂ / PI	ADULT / ADOLESCENT

DIMENSION	DISPLAY DIRECTION	POWER SUPPLY	WEIGHT
62 x 37 x 32 mm	4 DIRECTIONS	AAA x 2	42.5g

RANGE & ACCURACY			INDICATION
SpO ₂	70 - 100%	± 2 %	REAL-TIME BATTERY STATUS
PULSE RATE	25-250 BPM	± 2 BPM	VALUE LIMITED REMIND

Multiple Display Modes

Pulse Rate Measuring

SpO₂ Measuring

Brightness Adjustable

Alert Function

JPD-500B

DISPLAY	PARAMETERS	APPLICATION RANGE
HD LED	PR / SpO ₂	ADULT / ADOLESCENT

DIMENSION	POWER SUPPLY	WEIGHT
58 x 32 x 33.2 mm	AAA x 2	49.5g

RANGE & ACCURACY			INDICATION
SpO ₂	70 - 100%	± 2 %	LOW BATTERY REMIND
PULSE RATE	25-250 BPM	± 2 BPM	VALUE LIMITED REMIND

Pulse Rate Measuring

SpO₂ Measuring



JPD-500C

DISPLAY	PARAMETERS	APPLICATION RANGE
OLED	PR / SpO ₂	ADULT / ADOLESCENT

DIMENSION	DISPLAY DIRECTION	POWER SUPPLY	WEIGHT
58 x 32 x 33.2 mm	4 DIRECTIONS	AAA x 2	49.4g

RANGE & ACCURACY			INDICATION
SpO ₂	70 - 100%	± 2 %	REAL-TIME BATTERY STATUS
PULSE RATE	25-250 BPM	± 2 BPM	VALUE LIMITED REMIND

Multiple Display Modes

Pulse Rate Measuring

SpO₂ Measuring

Brightness Adjustable

Alert Function





JPD-500G

DISPLAY
HD LED

PARAMETERS
PR / SpO₂ / PI

APPLICATION RANGE
ADULT / ADOLESCENT

DIMENSION
58 x 32 x 32.9 mm

POWER SUPPLY
AAA x 2

WEIGHT
50.4g

RANGE & ACCURACY
SpO₂ 70 - 99% ± 2 %
PULSE RATE 25-250 BPM ± 2 BPM

INDICATION
LOW BATTERY REMIND
VALUE LIMITED REMIND



Pulse Rate
Measuring



SpO₂
Measuring



Brightness
Adjustable



Alert
Function



JPD-510E

DISPLAY
HD LED

PARAMETERS
PR / SpO₂

APPLICATION RANGE
ADULT / ADOLESCENT

DIMENSION
62 x 37 x 32 mm

POWER SUPPLY
AAA x 2

WEIGHT
50g

RANGE & ACCURACY
SpO₂ 70 - 100% ± 2 %
PULSE RATE 25-250 BPM ± 2 BPM

INDICATION
LOW BATTERY REMIND
VALUE LIMITED REMIND



Pulse Rate
Measuring



SpO₂
Measuring



Alert
Function



Rotating Display
(Optional)



JPD-500E

DISPLAY
OLED

PARAMETERS
PR / SpO₂ / PI

APPLICATION RANGE
ADULT / ADOLESCENT

DIMENSION
62 x 37 x 32 mm

DISPLAY DIRECTION
4 DIRECTIONS

POWER SUPPLY
AAA x 2

WEIGHT
42.5g

RANGE & ACCURACY
SpO₂ 70 - 100% ± 2 %
PULSE RATE 25-250 BPM ± 2 BPM

INDICATION
REAL-TIME BATTERY STATUS
VALUE LIMITED REMIND



Multiple
Display Modes



Pulse Rate
Measuring



SpO₂
Measuring



Brightness
Adjustable



Alert
Function



JPD-500H

DISPLAY
OLED

PARAMETERS
PR / SpO₂ / PI

APPLICATION RANGE
CHILD (>3 Years Old)

DIMENSION
52.6 x 34.1 x 31.2 mm

DISPLAY DIRECTION
4 DIRECTIONS (OLED)

POWER SUPPLY
AAA x 2

WEIGHT
48.4g

RANGE & ACCURACY
SpO₂ 70 - 100% ± 2 %
PULSE RATE 25-250 BPM ± 2 BPM

INDICATION
REAL-TIME BATTERY STATUS
VALUE LIMITED REMIND



Design
For Child



Pulse Rate
Measuring



SpO₂
Measuring



Brightness
Adjustable



Alert
Function



JPD-502F

DISPLAY
OLED

PARAMETERS
PR / SpO₂ / PI

APPLICATION RANGE
ADULT / ADOLESCENT

DIMENSION
62 x 37 x 32 mm

DISPLAY DIRECTION
4 DIRECTIONS

POWER SUPPLY
AAA x 2

WEIGHT
42.5g

RANGE & ACCURACY
SpO₂ 70 - 100% ± 2 %
PULSE RATE 25-250 BPM ± 2 BPM

INDICATION
REAL-TIME BATTERY STATUS
VALUE LIMITED REMIND



Multiple
Display Modes



Pulse Rate
Measuring



SpO₂
Measuring



Brightness
Adjustable



Alert
Function



JPD-510D

DISPLAY HD LED	PARAMETERS PR / SpO ₂	APPLICATION RANGE ADULT / ADOLESCENT
DIMENSION 62.2 x 37.0 x 33.1 mm	POWER SUPPLY AAA x 2	WEIGHT (WITHOUT BATTERIES) 28g
RANGE & ACCURACY SpO ₂ 70 - 100% ± 2 % PULSE RATE 25-250 BPM ± 2 BPM		INDICATION LOW BATTERY REMIND VALUE LIMITED REMIND


Pulse Rate Measuring


SpO₂ Measuring


Alert Function


Rotating Display (Optional)



JPD-500D

DISPLAY OLED	PARAMETERS PR / SpO ₂ / PI	APPLICATION RANGE ADULT / ADOLESCENT
DIMENSION 62.2 x 37.0 x 33.1 mm	DISPLAY DIRECTION 4 DIRECTIONS	POWER SUPPLY AAA x 2
WEIGHT 42.5g		
RANGE & ACCURACY SpO ₂ 70 - 100% ± 2 % PULSE RATE 25-250 BPM ± 2 BPM		INDICATION REAL-TIME BATTERY STATUS VALUE LIMITED REMIND


Multiple Display Modes


Pulse Rate Measuring


SpO₂ Measuring


Brightness Adjustable


Alert Function


Rotating Display (Optional)





JPD-500A

DISPLAY OLED	PARAMETERS PR / SpO ₂	APPLICATION RANGE ADULT / ADOLESCENT
DIMENSION 64 x 35 x 34 mm	DISPLAY DIRECTION 4 DIRECTIONS	POWER SUPPLY AAA x 2
WEIGHT 57g		
RANGE & ACCURACY SpO ₂ 70 - 100% ± 2 % PULSE RATE 25-250 BPM ± 2 BPM		INDICATION REAL-TIME BATTERY STATUS VALUE LIMITED REMIND


Multiple Display Modes


Pulse Rate Measuring


SpO₂ Measuring


Brightness Adjustable


Alert Function

Fingertip Pulse Oximeter	      					
	Parameters					
 Product Model	JPD-500G	SPA11	JPD-501F	JPD-500F	SPA10	SPA20
 Display Screen	LED	TFT	OLED		LED	
 PI	✓				—	
 Battery Status	Low Battery Remind	Real Time			Low Battery Remind	
 Plethysmograph Waveform	—	✓			—	
 Directions	1	4			1	2
 Display Modes	1	6			1	2
 Brightness Control	✓	—	✓		—	
 Bluetooth	✕				—	
 Rotating Display (Optional)	—		✓	—		
 Application Range	ADULT / ADOLESCENT					
 Dimension (mm)	58 x 32 x 32.9	60 x 33 x 32	62 x 37 x 32	62 x 37 x 32	62 x 33 x 32	67 x 35 x 31.7
 Weight(g) (Include Batteries)	50.4	49.5	42.5	42.5	49.5	60
 Functions	SpO ₂ PR PI Limit Remind	SpO ₂ PR PI Limit Remind Demo Mode	SpO ₂ PR PI Limit Adjust & Remind Demo Mode		SpO ₂ PR Limit Remind	

Fingertip Pulse Oximeter Parameters													
 Product Model	SPA30	JPD-500B	JPD-510D	JPD-510E	JPD-500H	JPD-500I		JPD-500A	JPD-500C	JPD-500D	JPD-500E	JPD-502F	JPD-502G
 Display Screen	LED				OLED	TFT		OLED					
 PI	—				✓			—		✓			
 Battery Status	Low Battery Remind				Real Time			Real Time					
 Plethysmograph Waveform	—				✓			✓					
 Directions	2	1			4			4					
 Display Modes	2	1			6			6					
 Brightness Control	—				✓	—		✓					
 Bluetooth	—							—					
 Rotating Display (Optional)	—				—			✓	—	✓	—	—	
 Application Range	ADULT / ADOLESCENT				CHILD (>3 YEARS OLD)	ADULT / ADOLESCENT		ADULT / ADOLESCENT					
 Dimensions (mm)	62 x 33 x 31.5	58 x 32 x 33.2	62.2 x 37.0 x 33.1	62 x 37 x 32	52.6 x 34.1 x 31.2	58 X 33.9 X 32		64 X 35 X 34	58 x 32 x 33.2	62.2 x 37.0 x 33.1	62 x 37 x 32	62 X 37 X 32	58.2 x 32.2 x 32
 Weight(g) (Include Batteries)	51.5	49.5	28 (without batteries)	50	48.4	49.5		57	49.4	42.5	42.5	42.5	51.5
 Functions	SpO ₂ PR Limit Remind				SpO ₂ PR PI Limit Adjust & Remind	SpO ₂ PR PI Limit Remind Demo Mode		SpO ₂ PR Demo Mode Limit Adjust & Remind		SpO ₂ PR PI Limit Adjust & Remind Demo Mode			

INFRARED THERMOMETER



Forehead/Ear
Auto Switch



Magnetic
Probe Cap



Child
Mode



10s Auto
Shut Off



Smart
Fever Alert

ONE KEY MEASUREMENT, ONE SECOND READING.

Infrared Thermometer uses an infrared sensor to detect the infrared radiation emitted by the eardrum and forehead, and converts the value through a microcomputer to obtain the temperature data.

Users can quickly get measurement results after positioning properly the temperature probe in the ear or on the forehead.

In addition to adults' daily temperature measurement, it can also effectively avoid measurement difficulties caused by crying resistance when measuring infants and young children. This kind of thermometer is more suitable for the whole family to use together.

FEATURES

- Child Mode
 - Silent Mode
 - Forehead /Ear/Object Mode
- Professional Accuracy
 - Hygienic Non-Touch
 - 1 Second Reading
- 10s Auto Shut Off
 - Smart Fever Alert
 - 20 Sets Records



NEW
ARRIVAL



JPD-FR302

DISPLAY HD LCD	MEMORY STORAGE 20 SETS DATA	TEMPERATURE UNITS °C/°F
AUTOMATIC SHUTDOWN 10S ± 1S	BATTERY AAA X 2	ACCURACY(CLINICAL) 0.3°C / 0.6°F
MEASURING RANGE		
FOREHEAD MODE	22 - 43 °C / 71.6 - 109.4 °F	
EAR MODE	34 - 43 °C / 93.2 - 109.4 °F	
OBJECT MODE	0 - 100 °C / 32.0 - 212.0 °F	
Forehead/Ear Auto Switch	Magnetic Probe Cap	Child Mode
10s Auto Shut Off	Smart Fever Alert	Silent Mode

NEW
ARRIVAL



JPD-FR304

DISPLAY HD LCD	MEMORY STORAGE 20 SETS DATA	TEMPERATURE UNITS °C/°F
AUTOMATIC SHUTDOWN 10S ± 1S	BATTERY AAA X 2	ACCURACY(CLINICAL) 0.3°C / 0.6°F
MEASURING RANGE		
EAR MODE	34 - 43 °C / 93.2 - 109.4 °F	
Magnetic Probe Cap	Child Mode	10s Auto Shut Off
Smart Fever Alert	Silent Mode	

NEW
ARRIVAL



JPD-FR403

DISPLAY HD LCD	MEMORY STORAGE 20 SETS DATA	TEMPERATURE UNITS °C/°F
AUTOMATIC SHUTDOWN 10S ± 1S	BATTERY AAA X 2	ACCURACY(CLINICAL) 0.3°C / 0.6°F
MEASURING RANGE		
FOREHEAD MODE	22 - 43 °C / 71.6 - 109.4 °F	
EAR MODE	34 - 43 °C / 93.2 - 109.4 °F	
OBJECT MODE	0 - 100 °C / 32.0 - 212.0 °F	
Forehead/Ear Auto Switch	Magnetic Probe Cap	Child Mode
10s Auto Shut Off	Smart Fever Alert	Silent Mode

NEW ARRIVAL



JPD-FR301

DISPLAY HD LED	MEMORY STORAGE 20 SETS DATA	TEMPERATURE UNITS °C/°F
AUTOMATIC SHUTDOWN 10S ± 1S	BATTERY AAA X 2	ACCURACY(CLINICAL) 0.3°C / 0.6°F
MEASURING RANGE		
FOREHEAD MODE	22 - 43 °C / 71.6 - 109.4 °F	
EAR MODE	34 - 43 °C / 93.2 - 109.4 °F	
OBJECT MODE	0 - 100 °C / 32.0 - 212.0 °F	
 Forehead/Ear Auto Switch	 Magnetic Probe Cap	 Child Mode
 10s Auto Shut Off	 Smart Fever Alert	 Silent Mode

NEW ARRIVAL



JPD-FR204

DISPLAY HD LED	MEMORY STORAGE 20 SETS DATA	TEMPERATURE UNITS °C/°F
AUTOMATIC SHUTDOWN 10S ± 1S	BATTERY AAA X 2	ACCURACY(CLINICAL) 0.3°C / 0.6°F
MEASURING RANGE		
FOREHEAD MODE	22 - 43 °C / 71.6 - 109.4 °F	
OBJECT MODE	0 - 100 °C / 32.0 - 212.0 °F	
 Hygienic Non-Touch	 Child Mode	 10s Auto Shut Off
 Smart Fever Alert	 Silent Mode	



JPD-FR410

DISPLAY HD LCD	MEMORY STORAGE 20 SETS DATA	TEMPERATURE UNITS °C/°F
AUTOMATIC SHUTDOWN 10S ± 1S	BATTERY AAA X 2	ACCURACY(CLINICAL) 0.3°C / 0.6°F
MEASURING RANGE		
FOREHEAD MODE	22 - 43 °C / 71.6 - 109.4 °F	
EAR MODE	34 - 43 °C / 93.2 - 109.4 °F	
OBJECT MODE	0 - 100 °C / 32.0 - 212.0 °F	
 Child Mode	 10s Auto Shut Off	 Smart Fever Alert
 Silent Mode		





JPD-FR409-BT

DISPLAY HD LCD	MEMORY STORAGE 20 SETS DATA	TEMPERATURE UNITS °C/°F
AUTOMATIC SHUTDOWN 10S ± 1S	BATTERY AAA X 2	ACCURACY(CLINICAL) 0.3°C / 0.6°F
MEASURING RANGE		
FOREHEAD MODE	22 - 43 °C / 71.6 - 109.4 °F	
EAR MODE	34 - 43 °C / 93.2 - 109.4 °F	
OBJECT MODE	0 - 100 °C / 32.0 - 212.0 °F	
 Magnetic Probe Cap	 Child Mode	 10s Auto Shut Off
 Smart Fever Alert	 Silent Mode	



JPD-FR409

DISPLAY HD LCD	MEMORY STORAGE 20 SETS DATA	TEMPERATURE UNITS °C/°F
AUTOMATIC SHUTDOWN 10S ± 1S	BATTERY AAA X 2	ACCURACY(CLINICAL) 0.3°C / 0.6°F
MEASURING RANGE		
FOREHEAD MODE	22 - 43 °C / 71.6 - 109.4 °F	
EAR MODE	34 - 43 °C / 93.2 - 109.4 °F	
OBJECT MODE	0 - 100 °C / 32.0 - 212.0 °F	
 Magnetic Probe Cap	 Child Mode	 10s Auto Shut Off
 Smart Fever Alert	 Silent Mode	





JPD-FR400

DISPLAY HD LCD	MEMORY STORAGE 20 SETS DATA	TEMPERATURE UNITS °C/°F
AUTOMATIC SHUTDOWN 10S ± 1S	BATTERY AAA X 2	ACCURACY(CLINICAL) 0.3°C / 0.6°F
MEASURING RANGE		
FOREHEAD/EAR MODE	32 - 42.2 °C / 89.6 - 107.9 °F	
OBJECT MODE	0.0 - 100 °C / 32.0 - 212.0 °F	
 1S	 10s Auto Shut Off	 Smart Fever Alert
1 Second Reading		



JPD-FR300

DISPLAY HD LCD	MEMORY STORAGE 20 SETS DATA	TEMPERATURE UNITS °C/°F
AUTOMATIC SHUTDOWN 10S ± 1S	BATTERY AAA X 2	ACCURACY(CLINICAL) 0.3°C / 0.6°F
MEASURING RANGE		
FOREHEAD MODE	22 - 43 °C / 71.6 - 109.4 °F	
EAR MODE	34 - 43 °C / 93.2 - 109.4 °F	
OBJECT MODE	0 - 100 °C / 32.0 - 212.0 °F	
 Forehead/Ear Auto Switch	 Magnetic Probe Cap	 Child Mode
 10s Auto Shut Off	 Smart Fever Alert	 Silent Mode



JPD-FR100+

DISPLAY HD LCD	MEMORY STORAGE 20 SETS DATA	TEMPERATURE UNITS °C/°F
AUTOMATIC SHUTDOWN 10S ± 1S	BATTERY AAA X 2	ACCURACY(CLINICAL) 0.3°C / 0.6°F
MEASURING RANGE		
FOREHEAD/EAR MODE	32 - 42.2 °C / 89.6 - 107.9 °F	
OBJECT MODE	0.0 - 100 °C / 32.0 - 212.0 °F	
 10s Auto Shut Off	 Smart Fever Alert	 1S Reading

NEW ARRIVAL



JPD-FR415

DISPLAY HD LCD	MEMORY STORAGE 20 SETS DATA	TEMPERATURE UNITS °C/°F
AUTOMATIC SHUTDOWN 10S ± 1S	BATTERY AAA X 2	ACCURACY(CLINICAL) 0.3°C / 0.6°F
MEASURING RANGE		
EAR MODE	34 - 43 °C / 93.2 - 109.4 °F	
 Magnetic Probe Cap	 Child Mode	 10s Auto Shut Off
 Smart Fever Alert	 Silent Mode	

NEW ARRIVAL



JPD-FR205

DISPLAY HD LCD	MEMORY STORAGE 20 SETS DATA	TEMPERATURE UNITS °C/°F
AUTOMATIC SHUTDOWN 10S ± 1S	BATTERY AAA X 2	ACCURACY(CLINICAL) 0.3°C / 0.6°F
MEASURING RANGE		
FOREHEAD MODE	22 - 43 °C / 71.6 - 109.4 °F	
OBJECT MODE	0 - 100 °C / 32.0 - 212.0 °F	
 Distance Detection	 Magnetic Battery Cap	 Hygienic Non-Touch
 10s Auto Shut Off	 Smart Fever Alert	 Silent Mode



JPD-FR401

DISPLAY HD LCD	MEMORY STORAGE 20 SETS DATA	TEMPERATURE UNITS °C/°F
AUTOMATIC SHUTDOWN 10S ± 1S	BATTERY AAA X 2	ACCURACY(CLINICAL) 0.3°C / 0.6°F
MEASURING RANGE		
FOREHEAD/EAR MODE	32 - 42.2 °C / 89.6 - 107.9 °F	
OBJECT MODE	0.0 - 100 °C / 32.0 - 212.0 °F	
 1S Reading	 10s Auto Shut Off	 Smart Fever Alert

NEW ARRIVAL



FTA10

DISPLAY HD LCD	MEMORY STORAGE 20 SETS DATA	TEMPERATURE UNITS °C/°F
AUTOMATIC SHUTDOWN 10S ± 1S	BATTERY AAA X 2	ACCURACY(CLINICAL) 0.3°C / 0.6°F
MEASURING RANGE		
FOREHEAD MODE	22 - 43 °C / 71.6 - 109.4 °F	
OBJECT MODE	0 - 100 °C / 32.0 - 212.0 °F	
 Hygienic Non-Touch	 Child Mode	 10s Auto Shut Off
 Smart Fever Alert	 Silent Mode	



JPD-FR203

DISPLAY HD LCD	MEMORY STORAGE 20 SETS DATA	TEMPERATURE UNITS °C/°F
AUTOMATIC SHUTDOWN 10S ± 1S	BATTERY AAA X 2	ACCURACY (CLINICAL) 0.3°C / 0.6°F
MEASURING RANGE		
FOREHEAD MODE	22 - 43 °C / 71.6 - 109.4 °F	
OBJECT MODE	0 - 100 °C / 32.0 - 212.0 °F	
Hygienic Non-Touch	Child Mode	10s Auto Shut Off
	Smart Fever Alert	Silent Mode



JPD-FR202

DISPLAY HD LCD	MEMORY STORAGE 20 SETS DATA	TEMPERATURE UNITS °C/°F
AUTOMATIC SHUTDOWN 10S ± 1S	BATTERY AAA X 2	ACCURACY (CLINICAL) 0.3°C / 0.6°F
MEASURING RANGE		
FOREHEAD MODE	22 - 43 °C / 71.6 - 109.4 °F	
OBJECT MODE	0 - 100 °C / 32.0 - 212.0 °F	
Hygienic Non-Touch	10s Auto Shut Off	Smart Fever Alert
	Silent Mode	

Infrared Thermometer Parameters												
Product Model		JPD-FR400		JPD-FR409-BT		JPD-FR403		JPD-FR302		JPD-FR301		
Display Screen		HD LCD									HD LED	
Forehead Mode		32.0-42.2 °C 89.6-107.9 °F		22.0-43.0 °C / 71.6-109.4 °F								
Ear Mode		0.0-100.0 °C 32.0-212.0 °F		34.0-43.0 °C / 93.2-109.4 °F								
Object Mode		0.0-100.0 °C / 32.0-212.0 °F										
Child Mode		—										
Silent Mode		—										
Backlight/ Indicator										—		
Clinical Accuracy		± 0.3 °C / ± 0.6 °F										
Auto Shut Off		10s ± 1s										
Features		❶ Fever warning ❷ °C/°F adjustable ❸ Memory storage: 20 sets data ❹ Batteries: AAA × 2 ❺ Forehead mode: for all ages ❻ Ear mode for children: ≥ 3 months old ❼ Child mode: ≤ 12 years old										

BLOOD PRESSURE MONITOR



- Bluetooth
- Guest Mode
- 2 Users
- Voice Switch
- Type-C
- 199 Records
- Average Value
- 30s Auto Power Off
- Irregular Heartbeat Indicator

FIND THE RIGHT BLOOD PRESSURE MONITOR FOR YOU.

Our blood pressure monitor features advanced technology to automatically inflate and deflate at the appropriate level, delivering a personalized accurate, and comfortable measurement.

It can also calculate average blood pressure value based on the last three readings. It is convenient for you to use at home or take it out because of the compact size, so you’ ll never miss a measurement.

FEATURES

- Universal cuff
- Extra-large screen
- One touch operation
- Display date and time
- Lightweight & portable
- Average of last 3 readings
- Irregular heartbeat indicator
- Auto-off power saving feature
- 2-User mode & memory function



NEW ARRIVAL



JPD-HA201C

TYPE	LCD SIZE	CUFF LENGTH(cm)
ARM TYPE	59 x 47 mm / 3 inch	22-36 / 22-42
USER MEMORY	AUTOMATIC SHUTDOWN	
2 x 199 RECORDS	30s	
MEASUREMENT RANGE	POWER SUPPLY	
PRESSURE VALUE	0-295 ± 3 mmHg	AAA x 4
PULSE RATE	40-199 ± 5% BPM	USB 5V
Average Value	Type-C	Display Date & Time
		2 Users
		Backlight (Optional)
		Irregular Heartbeat Indicator

NEW ARRIVAL

AFib



JPD-HA200C

TYPE	LCD SIZE	CUFF LENGTH(cm)
ARM TYPE	72 x 58 mm / 3.6 inch	22-36 / 22-42
USER MEMORY	AUTOMATIC SHUTDOWN	
2 x 199 RECORDS	30s	
MEASUREMENT RANGE	POWER SUPPLY	
PRESSURE VALUE	0-295 ± 3 mmHg	AAA x 4
PULSE RATE	40-199 ± 5% BPM	USB 5V
Guest Mode	Average Value	Type-C
		2 Users
		Backlight (Optional)
		Irregular Heartbeat Indicator

Bluetooth



HAA11

BLUETOOTH	WITHOUT BLUETOOTH	TYPE	LCD SIZE
HAA11	HAA10	ARM TYPE	91 x 77 mm / 4.8 inch
USER MEMORY	AUTOMATIC SHUTDOWN	CUFF LENGTH(cm)	
2 x 199 RECORDS	30s	22-36 / 22-42	
MEASUREMENT RANGE	POWER SUPPLY		
PRESSURE VALUE	0-295 ± 3 mmHg	AA x 4	
PULSE RATE	40-199 ± 5% BPM	USB 5V	
Guest Mode	Voice Switch	Type-C	2 Users
		Average Value	Irregular Heartbeat Indicator



HAA21

BLUETOOTH HAA21	WITHOUT BLUETOOTH HAA20	TYPE ARM TYPE	LCD SIZE 71 x 101 mm / 4.9 inch
--------------------	----------------------------	------------------	------------------------------------

USER MEMORY 2 x 199 RECORDS	AUTOMATIC SHUTDOWN 30s	CUFF LENGTH(cm) 22-36 / 22-42
--------------------------------	---------------------------	----------------------------------

MEASUREMENT RANGE PRESSURE VALUE PULSE RATE	0-295 ± 3 mmHg 40-199 ± 5% BPM	POWER SUPPLY AAA x 4 USB 5V
---	-----------------------------------	-----------------------------------



Guest Mode



Voice Switch



Type-C



Average Value



2 Users



Irregular Heartbeat Indicator



JPD-HA121

TYPE ARM TYPE	LCD SIZE 83 x 97 mm / 5.1 inch	CUFF LENGTH(cm) 22-36 / 22-42
------------------	-----------------------------------	----------------------------------

USER MEMORY 2 x 99 RECORDS	AUTOMATIC SHUTDOWN 30s
-------------------------------	---------------------------

MEASUREMENT RANGE PRESSURE VALUE PULSE RATE	0-295 ± 3 mmHg 40-199 ± 5% BPM	POWER SUPPLY AA x 4 USB 5V
---	-----------------------------------	----------------------------------



Extra-Large Screen



Average Value



30s Auto Power Off



2 Users



Backlight



Irregular Heartbeat Indicator



JPD-HA120

TYPE ARM TYPE	LCD SIZE 60 x 77 mm / 3.9 inch	CUFF LENGTH(cm) 22-36 / 22-42
------------------	-----------------------------------	----------------------------------

USER MEMORY 2 x 99 RECORDS	AUTOMATIC SHUTDOWN 30s
-------------------------------	---------------------------

MEASUREMENT RANGE PRESSURE VALUE PULSE RATE	0-295 ± 3 mmHg 40-199 ± 5% BPM	POWER SUPPLY AA x 4 USB 5V
---	-----------------------------------	----------------------------------



Average Value



30s Auto Power Off



Display Date & Time



2 Users



Irregular Heartbeat Indicator



JPD-HA100

TYPE ARM TYPE	LCD SIZE 60 x 77 mm / 3.9 inch	CUFF LENGTH(cm) 22-36 / 22-42
------------------	-----------------------------------	----------------------------------

USER MEMORY 2 x 99 RECORDS	AUTOMATIC SHUTDOWN 30s
-------------------------------	---------------------------

MEASUREMENT RANGE PRESSURE VALUE PULSE RATE	0-295 ± 3 mmHg 40-199 ± 5% BPM	POWER SUPPLY AA x 4 USB 5V
---	-----------------------------------	----------------------------------



Average Value



30s Auto Power Off



Display Date & Time



2 Users



Irregular Heartbeat Indicator

B

JPD-HA100



TYPE	LCD SIZE	CUFF LENGTH(cm)
ARM TYPE	60 x 77 mm / 3.9 inch	22-36 / 22-42
USER MEMORY	AUTOMATIC SHUTDOWN	
2 x 99 RECORDS	30s	
MEASUREMENT RANGE		POWER SUPPLY
PRESSURE VALUE	0-295 ± 3 mmHg	AA x 4
PULSE RATE	40-199 ± 5% BPM	USB 5V

- Average Value
- 30s Auto Power Off
- Display Date & Time
- 2 Users
- Backlight
- Irregular Heartbeat Indicator

JPD-HA210



TYPE	LCD SIZE	CUFF LENGTH(cm)
ARM TYPE	59 x 47 mm / 3 inch	22-36 / 22-42
USER MEMORY	AUTOMATIC SHUTDOWN	
2 x 120 RECORDS	30s	
MEASUREMENT RANGE		POWER SUPPLY
PRESSURE VALUE	0-295 ± 3 mmHg	AAA x 4
PULSE RATE	40-199 ± 5% BPM	USB 5V

- Average Value
- Type-C
- Display Date & Time
- 2 Users
- Backlight (Optional)
- Irregular Heartbeat Indicator

JPD-HA300



TYPE	LCD SIZE	CUFF LENGTH(cm)
ARM TYPE	59 x 47 mm / 3 inch	22-36 / 22-42
USER MEMORY	AUTOMATIC SHUTDOWN	
2 x 120 RECORDS	30s	
MEASUREMENT RANGE		POWER SUPPLY
PRESSURE VALUE	0-295 ± 3 mmHg	AAA x 4
PULSE RATE	40-199 ± 5% BPM	USB 5V

- Average Value
- Type-C
- Display Date & Time
- 2 Users
- Three-color backlight
- Irregular Heartbeat Indicator

HWA20



TYPE	LCD SIZE	CUFF LENGTH(cm)
WRIST	41x 41 mm / 2.3 inch	13.5-21.5
USER MEMORY	POWER SUPPLY	AUTO SHUTDOWN
2 x 199 RECORDS	AAA x 2	30s
MEASUREMENT RANGE		
PRESSURE VALUE	0-295 ± 3 mmHg	
PULSE RATE	40-199 ± 5% BPM	

- Guest Mode
- Voice Switch
- Average Value
- 2 Users
- Irregular Heartbeat Indicator
- Display Date & Time

JPD-HA200



TYPE	LCD SIZE	CUFF LENGTH(cm)
ARM TYPE	59 x 47 mm / 3 inch	22-36 / 22-42
USER MEMORY	AUTOMATIC SHUTDOWN	
2 x 120 RECORDS	30s	
MEASUREMENT RANGE		POWER SUPPLY
PRESSURE VALUE	0-295 ± 3 mmHg	AAA x 4
PULSE RATE	40-199 ± 5% BPM	USB 5V

- Average Value
- Type-C
- Display Date & Time
- 2 Users
- Backlight (Optional)
- Irregular Heartbeat Indicator

HWA10



TYPE	LCD SIZE	CUFF LENGTH(cm)
WRIST	45 x 41 mm / 2.4 inch	13.5-21.5
USER MEMORY	POWER SUPPLY	AUTO SHUTDOWN
2 x 199 RECORDS	AAA x 2	30s
MEASUREMENT RANGE		
PRESSURE VALUE	0-295 ± 3 mmHg	
PULSE RATE	40-199 ± 5% BPM	

- Guest Mode
- Voice Switch
- Average Value
- 2 Users
- Irregular Heartbeat Indicator
- Display Date & Time





HWA11

TYPE	LCD SIZE	CUFF LENGTH(cm)
WRIST	45 x 41 mm / 2.4 inch	13.5-21.5

USER MEMORY	POWER SUPPLY	AUTO SHUTDOWN
2 x 199 RECORDS	AAA x 2	30s

MEASUREMENT RANGE

PRESSURE VALUE

PULSE RATE

0-295 ± 3 mmHg

40-199 ± 5% BPM













Guest Mode

Voice Switch

Average Value

2 Users

Irregular Heartbeat Indicator

Display Date & Time





HWA21

TYPE	LCD SIZE	CUFF LENGTH(cm)
WRIST	41x 41 mm / 2.3 inch	13.5-21.5

USER MEMORY	POWER SUPPLY	AUTO SHUTDOWN
2 x 199 RECORDS	AAA x 2	30s

MEASUREMENT RANGE

PRESSURE VALUE

PULSE RATE

0-295 ± 3 mmHg

40-199 ± 5% BPM













Guest Mode

Voice Switch


Average Value

2 Users

Irregular Heartbeat Indicator

Display Date & Time

Blood Pressure Monitor Parameters		    				
 Product Model	JPD-HA201C	JPD-HA200C	HAA11	HAA21	HAA10	
 Type	Arm					
 Pressure Value	0-295 ± 3 mmHg					
 Pulse Rate	40-199 ± 5% BPM					
 Display Size mm inch	59 x 47	72 x 58	91 x 77	71 x 101	91 x 77	
 Average						
 Bluetooth						
 Backlight	1 2 3 4  —					
 Auto Power Off	30s					
 Cuff Length (cm)	22-36 / 22-42 (Optional)					
 Power Supply	AAA x 4 USB 5V		AA x 4 USB 5V	AAA x 4 USB 5V	AA x 4 USB 5V	
 User Memory	2 x 199 records					
 Type-C						

Blood Pressure Monitor Parameters													
 Product Model	HAA20	JPD-HA300	JPD-HA210	JPD-HA200	JPD-HA120			JPD-HA121	JPD-HA100	JPD-HA100	JPD-HA101(LCD)	JPD-HA101(LED)	
 Type	Arm							Arm					
 Pressure Value	0-295 ± 3 mmHg							0-295 ± 3 mmHg					
 Pulse Rate	40-199 ± 5% BPM							40-199 ± 5% BPM					
 Display Size	mm	71 x 101	59 x 47		59 x 47	60 x 77			83 x 97	60 x 77		83 x 97	83 x 97
	inch	4.9	3		3	3.9			5.1	3.9		5.1	5.1
 Average													
 Bluetooth	—								—				
 Backlight						—				—			—
 Auto Power Off	30s							30s					
 Cuff Length (cm)	22-36 / 22-42 (Optional)							22-36 / 22-42 (Optional)					
 Power Supply	AAA x 4 USB 5V				AA x 4 USB 5V			AA x 4 USB 5V					
 User Memory	2 x 199 records	2 x 120 records			2 x 99 records			2 x 99 records					
 Type-C					—			—					

Blood Pressure Monitor Parameters					
Product Model	HWA11	HWA21	HWA10	HWA20	
Type	Wrist				
Pressure Value	0-295 ± 3 mmHg				
Pulse Rate	40-199 ± 5% BPM				
Display Size	mm 45 x 41 inch 2.4	41 x 41 2.3	45 x 41 2.4	41 x 41 2.3	
Average	✓				
Bluetooth	✕		—		
Backlight	☀				
Auto Power Off	30s				
Cuff Length (cm)	13.5-21.5				
Power Supply	AAA x 2				
User Memory	2 x 199 records				
Type-C	—				

TENS THERAPY DEVICE



- 
- 
APP
Control
- 
5 Modes
15 Levels
- 
One Button
Control
- 
Automatic
Shutdown
- 
Relieve Pain
- 
Rechargeable
Battery


RELIEVE PAIN, SOOTHE MIND.

Small, portable and powerful TENS devices from JUMPER offer drug-free pain relief for sore or aching muscles of the lower back, shoulder, arm, leg and more.

The unit also features one or two independent stimulation channels and four self-adhesive electrodes. It offers a wide range of functions for increasing general well-being, pain relief, maintaining physical fitness, relaxation, muscle revitalization and combating tiredness. For these purposes, you can either choose from pre-set programs or specify your own to suit your individual needs.

FEATURES

- Pain relief
 - Pad off Alert
 - Automatic shutdown
- Low battery indication
 - Ultra-thin & lightweight
 - Easy to charge and use
- Various massage methods
 - Powerful & effective program







TREATMENT TIME	PULSE FREQUENCY	INTENSITY
1-60 mins	1-120HZ	0-15
MAX OUTPUT INTENSITY	POWER SUPPLY TYPE	PULSE WIDTH
30V at 500Ω	LITHIUM	0.15-0.20ms
DIMENSION(mm)	MAIN UNIT WEIGHT	
MAIN UNIT	44 x 44 x 13	
ELECTRODE PAD	φ90 / φ60	



APP ControlOne Button ControlPad Off Alert5 ModesRechargeable BatteryLow Battery Indication



TREATMENT TIME	PULSE FREQUENCY	INTENSITY
5-60 mins	1-120HZ	0-15
MAX OUTPUT INTENSITY	POWER SUPPLY TYPE	PULSE WIDTH
48mA at 500Ω	DRY BATTERIES	0.15-0.20ms
DIMENSION(mm)	MAIN UNIT WEIGHT	
MAIN UNIT	121.8 x 63.6 x 22.2	
ELECTRODE PAD	50 x 50	



5 Modes15 Levels IntensityPad Off AlertAvailable For 2 UsersLow Battery Indication





TREATMENT TIME	PULSE FREQUENCY	INTENSITY
5-60 mins	1-100HZ ± 10%	0-20
MAX OUTPUT INTENSITY	POWER SUPPLY TYPE	PULSE WIDTH
50V at 500Ω	LITHIUM	80μs
DIMENSION(mm)	MAIN UNIT WEIGHT	
MAIN UNIT	138x 68.7 x 13	
ELECTRODE PAD	50 x 50	




Type-CRechargeable Battery20 Modes20 Levels IntensityAvailable For 2 UsersPad Off Alert



A


JPD-ES220

TREATMENT TIME	PULSE FREQUENCY	INTENSITY
5-60 mins	1-100HZ ± 10%	0-20
MAX OUTPUT INTENSITY	POWER SUPPLY TYPE	PULSE WIDTH
50V at 500Ω	LITHIUM	80μs
DIMENSION(mm)	MAIN UNIT WEIGHT	
MAIN UNIT	138x 68.7 x 13	
ELECTRODE PAD	50 x 50	



JPD-ES230


TREATMENT TIME	PULSE FREQUENCY	INTENSITY
5-60 mins	2-120HZ	0-40
MAX OUTPUT INTENSITY	POWER SUPPLY TYPE	PULSE WIDTH
50V at 500Ω	DRY BATTERIES	0.05-0.30ms
DIMENSION(mm)	MAIN UNIT WEIGHT	
MAIN UNIT	130 x 64 x 23	
ELECTRODE PAD	50 x 50	
98g (WITHOUT BATTERY)		



A

JPD-ES210

TREATMENT TIME	PULSE FREQUENCY	INTENSITY
5-60 mins	1-120HZ	0-15
MAX OUTPUT INTENSITY	POWER SUPPLY TYPE	PULSE WIDTH
48mA at 500Ω	DRY BATTERIES	0.1-0.2ms
DIMENSION(mm)	MAIN UNIT WEIGHT	
MAIN UNIT	132.3 x 64.8 x 19.6	
ELECTRODE PAD	98 x 64	
	105g (WITHOUT BATTERY)	



LPA11


TREATMENT TIME	PULSE FREQUENCY	INTENSITY
5-60 mins	2-100HZ	0-40
MAX OUTPUT INTENSITY	POWER SUPPLY TYPE	PULSE WIDTH
50V at 500Ω	DRY BATTERIES	0.05-0.30ms
DIMENSION(mm)	MAIN UNIT WEIGHT	
MAIN UNIT	113x 62 x 22.5	
ELECTRODE PAD	50 x 50	



B


JPD-ES210

TREATMENT TIME	PULSE FREQUENCY	INTENSITY
5-60 mins	2-120HZ	0-15
MAX OUTPUT INTENSITY	POWER SUPPLY TYPE	PULSE WIDTH
48mA at 500Ω	DRY BATTERIES	0.05-0.30ms
DIMENSION(mm)	MAIN UNIT WEIGHT	
MAIN UNIT	132.3 x 64.8 x 19.6	
ELECTRODE PAD	98 x 64	
	105g (WITHOUT BATTERY)	



LPA12

TREATMENT TIME	PULSE FREQUENCY	INTENSITY
5-60 mins	2-120HZ	0-50
MAX OUTPUT INTENSITY	POWER SUPPLY TYPE	PULSE WIDTH
50V at 500Ω	LITHIUM	0.05-0.30ms
DIMENSION(mm)	MAIN UNIT WEIGHT	
MAIN UNIT	151 x 78 x 32	
ELECTRODE PAD	50 x 50	
	147g (WITHOUT BATTERY)	



★NEW★
ARRIVAL

TENS

MASSAGE

EMS

FITNESS



LPB12

TREATMENT TIME	PULSE FREQUENCY	INTENSITY
5-60 mins	2-120HZ	0-50
MAX OUTPUT INTENSITY	POWER SUPPLY TYPE	PULSE WIDTH
50V at 500Ω	LITHIUM	0.05-0.30ms
DIMENSION(mm)	MAIN UNIT WEIGHT	
MAIN UNIT	128g (WITHOUT BATTERY)	
ELECTRODE PAD	50 x 50	

79 Programs

50 Levels Intensity

Type-C

Pad Off Alert

Available For 2 Users

Low Battery Indication

TENS Therapy Device Parameters	 B  A  A   Bluetooth				
 Product Model	JPD-ES220	JPD-ES220	JPD-ES210	JPD-ES200	JPD-ES100
 TENS EMS MESSAGE	TENS - EMS - MESSAGE	TENS			
 Treatment Time	5-60 mins		5-60 mins	5-60 mins	1-60 mins
 Power Supply Type	Lithium		Dry batteries	Dry batteries	Lithium
 Max Output Intensity (at 500Ω)	50 V		48 mA	48 mA	30 V
 Pulse Width	80μs		0.1-0.2ms	0.15-0.20ms	
 Pulse Frequency	1-100Hz		1-120Hz		
 Low Battery Indication	✓				
 Number of Channels	2				1
 Auto Power Off	✓				
 Bluetooth	—				
 Main Unit Size (mm)	138 x 68.7 x 13		132.3 x 64.8 x 19.6	121.8 x 63.6 x 22.2	44 x 44 x 13
 Electrode Pad Size (mm)	50 x 50		98 x 64	50 x 50	φ90 / φ60
 Main Unit Weight (g)	123.9		105	86	50
 Modes	25	20	14	5	

PELVIC FLOOR REHABILITATION DEVICE



-  Bluetooth
-  HD LCD
-  2 separate channels
-  Type-C Charging
-  18 treatment modes
-  Customized mode

PORTABLE POSTPARTUM REHABILITATION THERAPY DEVICE

BFC10 is a portable Postpartum Rehabilitation Therapy Device, specifically designed for postnatal and menopausal women. It is based on the theory of low frequency electric therapy and pressure biofeedback training, offering a variety of treatments that have been optimized for maximum effectiveness. By stimulate the specific positions of the human body, BFC10 can help with urinary incontinence, repair rectus abdominis muscle, strengthen the vaginal and pelvic floor muscle and provide pelvic floor training to improve the health of women.



NEW
ARRIVAL



BFD10

OUTPUT CHANNELS	TREATMENT TIME	POWER SUPPLY
2	5-60 mins	Lithium Battery

NET WEIGHT	DIMENSION
145g	151mm×78mm×32mm

WORKING MODE
Continuous
operating

-  HD LCD
-  2 separate channels
-  Type-C Charging
-  18 treatment modes
-  Customized mode

NEW
ARRIVAL



BFC10

OUTPUT CHANNELS	TREATMENT TIME	POWER SUPPLY
2	5-60 mins	Lithium Battery

NET WEIGHT	DIMENSION
128g	151mm×82mm×33mm

WORKING MODE
Continuous
operating

-  HD LCD
-  2 separate channels
-  Type-C Charging
-  18 treatment modes
-  Customized mode

Pelvic Floor Rehabilitation Device Parameters			
 Product Model	BFD10	BFC10	
 Output Channels	2		
 Treatment Time	5-60 mins		
 Power Supply	Lithium Battery		
 Auto Power Off	✓		
 Low Power Indication	✓		
 Dimension (mm)	151 x 78 x 32	151 x 82 x 33	
 Net Weight (g)	145	128	
 Accessories	Vaginal Electrode Vaginal Pressure Biofeedback Probe Breast Electrode Electrode (100mmx50mm) Electrode (50mmX50mm)	Vaginal Electrode Vaginal Pressure Biofeedback Probe Breast Electrode Electrode (100mmx50mm) Electrode (50mmX50mm)	
 Modes	18		
 Working Mode	Continuous operating		

MESH NEBULIZER

COMPRESSOR

NEBULIZER



- 
Lightweight
and Portable
- 
Well
absorbed
- 
High
efficiency
- 
Low noise
- 
Independent
nebulizing
components

COMFORTABLE NEBULIZATION

Our nebulizer is high-performance and convenient. It fits various types of nebulized medication and can nebulize large particles of medication and viscous liquids. The nebulizer offers an affordable and efficient alternative for patients experiencing asthma, COPD or other respiratory disease conditions.



MTA10

POWER	MEDICATION CAPACITY	NEBULIZATION RATE
DC 3V	14 ML MAXIMUM	≥ 0.2 ML/MIN

SOUND LEVEL	PARTICLE SIZE
≤ 45 DB	1-5µm≥ 50 % (MMAD:3.5µm)

DIMENSION(mm)	Net Weight
MAIN UNIT	104g

- 
Lightweight
and Portable
- 
One-button
Operation
- 
Type-C / Battery
Power Supply
- 
Low noise
- 
14 ml Large
Capacity
- 
Self-Cleaning



CTA10

POWER	MEDICATION CAPACITY	NEBULIZATION RATE
AC 220V	8 ML MAXIMUM	≥ 0.2 ML/MIN

SOUND LEVEL	PARTICLE SIZE
≤ 50 DB	3.5 µm ±25%

DIMENSION(mm)	Net Weight
MAIN UNIT	1450g

- 
Easy to use
- 
Well
absorbed
- 
Adjustable
atomization
- 
Low noise
- 
Effective
and fast
medication
delivery
- 
More fully



CTB10

POWER	MEDICATION CAPACITY	NEBULIZATION RATE
AC 220V	8 ML MAXIMUM	≥ 0.2 ML/MIN

SOUND LEVEL	PARTICLE SIZE
≤ 50 DB	3.0 µm ±10%

DIMENSION(mm)	Net Weight
MAIN UNIT	170g

- 
Lightweight
and portable
- 
Well
absorbed
- 
Adjustable
atomization
- 
Low noise
- 
Effective
and fast
medication
delivery
- 
More fully

Nebulizer Comparison Table Parameters				
Product Model	MTA10	CTA10	CTB10	
Product Type	Mesh Nebulizer	Compressor Nebulizer	Compressor Nebulizer	
Nebulization Rate	≥ 0.2ml/min	≥ 0.2ml/min	≥ 0.2ml/min	
Power Supply	DC 5V / AAx2	AC 220V/50Hz	AC 100-240V-50 / 60Hz 0.4A	
Particles	1-5μm ≥ 50% (MMAD: 3.5μm)	3.5μm±25%	3.0μm±10%	
Sound	≤ 45dB	≤ 50 dB	≤ 50 dB	
Dimension (mm)	62.1 x 45 x 125.5	167 x 167 x 102	100 x 100 x 49.5	
Net Weight (g)	104	1450	170	
Accessories	Nebulization cup Adult mask Child mask Mouthpiece Headstrap	Nebulization cup Adult mask Child mask Mouthpiece Headstrap Airtube Air filters	Nebulization cup Adult mask Child mask Mouthpiece Headstrap Airtube Air filters	
Applicable Users	Elderly Adults Children	Elderly Adults Children	Elderly Adults Children	
Applicable Scenarios	Family	Hospital Clinic Family	Hospital Clinic Family	

WEIGHT SCALE

BODY FAT SCALE

- 

BMI
- 

Body Fat Rate
- 

Body Fat
- 

Body Water
- 

Automatic Power-Off

PROFESSIONAL MEASUREMENT, SCIENTIFICALLY MANAGE HEALTH.

These products are designed to measure body weight and body fat rate.

It can be connected to a smart phone via Bluetooth, and the measured data will be transmitted to the phone synchronously, which means that users can check their weight anytime by turning on the APP.

These products help users to develop good habits of diet management and exercise for the purpose of body management.

FEATURES


- Extra-large platform.

•Elegant and stylish design.

•Extra-high weight capacity.
- Accurate weight measurement.

•Durable tempered glass surface.

•Large easy-to-read LED/LCD display.
- The data of measurement can be checked on APP (JUMPER-HEALTH) via Bluetooth.





JPD-BFS710

DISPLAY	DIMENSIONS	POWER SUPPLY	
HD LED	325 x 325 x 23.5 mm	AAA x 4	

NET WEIGHT	MAX CAPACITY	DIVISION	UNIT
3 KG	180 KG / 396 LB	0.1 KG / 0.2 LB	KG / LB / ST

FUNCTION

LOW POWER & OVERLOAD INDICATION

AUTO POWER OFF IN 10 SECONDS



ITO



BMI



Body Fat Rate



Body Fat



Body Water



Automatic Power-Off



JPD-BFS100

DISPLAY	DIMENSIONS	POWER SUPPLY	UNIT
HD LED	372 x 298 x 29 mm	AAA x 3	KG / LB

NET WEIGHT	MAX CAPACITY	DIVISION	
3 KG	180 KG / 396 LB	0.1 KG / 0.2 LB	

FUNCTION

LOW POWER & OVERLOAD INDICATION

AUTO POWER OFF IN 10 SECONDS

BMI



Body Fat Rate



Body Fat



Body Water



Automatic Power-Off





JPD-BFS200B



JPD-BFS200D

JPD-BFS200A

DISPLAY	DIMENSIONS	POWER SUPPLY	
HD LED	280 x 280 x 23.5 mm	AAA x 3	

NET WEIGHT	MAX CAPACITY	DIVISION	UNIT
1.3 KG	200 KG / 440 LB	0.1 KG / 0.2 LB	KG / LB

FUNCTION

LOW POWER & OVERLOAD INDICATION

AUTO POWER OFF IN 10 SECONDS

BMI



Body Fat Rate



Body Fat



Body Water



Automatic Power-Off



JPD-BFS300



DISPLAY	DIMENSIONS	UNIT	POWER SUPPLY
HD LCD	280 x 280 x 23.5 mm	KG / LB / ST	AAA x 4

NET WEIGHT	MAX CAPACITY	DIVISION
1.3 KG	200 KG / 440 LB	0.1 KG / 0.2 LB

FUNCTION
LOW POWER & OVERLOAD INDICATION
AUTO POWER OFF IN 10 SECONDS

BMI

Body Fat Rate

Body Fat

Body Water

Automatic Power-Off



JPD-700A



DISPLAY	DIMENSIONS	MAX CAPACITY
HD LED	375 x 300 x 24.5 mm	200 KG / 440 LB

NET WEIGHT	DIVISION	UNIT	POWER SUPPLY
2.1 KG	0.1 KG / 0.2 LB	KG / LB	AAA x 3

FUNCTION
LOW POWER & OVERLOAD INDICATION
AUTO POWER OFF IN 10 SECONDS

Low Battery Indication

Unit Conversion

Automatic Power-Off

LED Screen Display

Top Quality Material

A



JPD-BS300



DISPLAY	DIMENSIONS	MAX CAPACITY
HD LCD	280 x 280 x 23.5 mm	200 KG / 440 LB / 31 ST

NET WEIGHT	DIVISION	UNIT	POWER SUPPLY
1.3 KG	0.1 KG / 0.2 LB	KG / LB / ST	AAA x 3

FUNCTION
LOW POWER & OVERLOAD INDICATION
AUTO POWER OFF IN 10 SECONDS

Low Battery Indication

Unit Conversion

Automatic Power-Off

LED Screen Display

Top Quality Material



JPD-700A



DISPLAY	DIMENSIONS	MAX CAPACITY
HD LED	375 x 300 x 24.5 mm	200 KG / 440 LB

NET WEIGHT	DIVISION	UNIT
2.1 KG	0.1 KG / 0.2 LB	KG / LB

FUNCTION
LOW POWER & OVERLOAD INDICATION
AUTO POWER OFF IN 10 SECONDS

Low Battery Indication

Unit Conversion

Automatic Power-Off

LED Screen Display

Top Quality Material

B



JPD-BS200



JPD-BS201

DISPLAY	DIMENSIONS	MAX CAPACITY
HD LED	280 x 280 x 23.5 mm	200 KG / 440 LB

NET WEIGHT	DIVISION	UNIT	POWER SUPPLY
1.3 KG	0.1 KG / 0.2 LB	KG / LB	AAA x 3

FUNCTION
LOW POWER & OVERLOAD INDICATION
AUTO POWER OFF IN 10 SECONDS

LED Screen Display

Top Quality Material

Automatic Power-Off

JPD-700C



DISPLAY	DIMENSIONS	MAX CAPACITY
HD LED	315 x 283 x 24.5 mm	180 KG / 396 LB

NET WEIGHT	DIVISION	UNIT	POWER SUPPLY
1.6 KG	0.1 KG / 0.2 LB	KG / LB	AAA x 4

FUNCTION
LOW POWER & OVERLOAD INDICATION
AUTO POWER OFF IN 10 SECONDS


Low Battery Indication

Unit Conversion

LED Screen Display

Top Quality Material

Automatic Power-Off

Weight & Body Fat Scale Parameters								 A B				
	Product Model	JPD-BFS710	JPD-BFS100	JPD-BFS200A/B/D	JPD-BFS300	JPD-BS300	JPD-BS200/201		JPD-700A	JPD-700A	JPD-700C	
	Display	HD LED		HD LED	HD LCD	HD LCD	HD LED		HD LED			
	Bluetooth	✕							✕		—	
	Power Supply	AAA x 4	AAA x 3	AAA x 3	AAA x 4	AAA x 3			AAA x 3		AAA x 4	
Dimensions (mm)	325 x 325 x 23.5	372 x 298 x 29	280 x 280 x 23.5					375x300x24.5		315 x 283 x 24.5		
Net Weight (Kg)	3	3	1.3					2.4		1.6		
Auto Power Off	10s							10s				
Max Capacity	180 KG / 396 LB		200 KG / 440 LB					200 KG / 440 LB		180 KG / 396 LB		
Indication	LOW POWER & OVERLOAD							LOW POWER & OVERLOAD				
Back Side												

REMOTE PATIENT MONITORING DEVICES



CELLULAR BLOOD PRESSURE MONITOR

4G
NETWORK



JPD-HA100

TYPE	LCD SIZE	CUFF LENGTH(cm)
ARM TYPE	60 x 77 mm / 3.9 inch	22-36 / 22-42
AUTOMATIC SHUTDOWN	DIMENSION (MM)	WEIGHT
30s	138x120x59	About 285g
MEASUREMENT RANGE		POWER SUPPLY
PRESSURE VALUE	0-295 ± 3 mmHg	AA x 4
PULSE RATE	40-199 ± 5% BPM	USB 5V

FEATURES

- Clinically validated accuracy, suitable for household and hospital
- 4G+2G Dual-mode communication, remote care anytime, anywhere
- Encrypted date transmission, sync data automatically
- One-button measurement, no additional operations
- No need to download apps and pair devices



CELLULAR FINGERTIP PULSE OXIMETER

4G
NETWORK



JPD-500D

DISPLAY OLED	PARAMETERS PR / SpO ₂ / PI	APPLICATION RANGE ADULT / ADOLESCENT
DIMENSION 62.2 x 37.0 x 33.1 mm	POWER SUPPLY AAA x 2	WEIGHT 42.5g
RANGE & ACCURACY SpO ₂ 70 - 100% ± 2 % PULSE RATE 25-250 BPM ± 2 BPM		INDICATION REAL-TIME BATTERY STATUS VALUE LIMITED REMIND

FEATURES

- One button measurement, no additional action required Low power consumption
- Upgrade sensor technology for ultimate accuracy
- Auto-sync to the cloud with no more operation
- Professional medical device suitable for clinical monitoring



CELLULAR WEIGHT SCALE

4G
NETWORK



JPD-700A

DISPLAY HD LCD	DIMENSIONS 375 x 300 x 29 mm	MAX CAPACITY 200 KG / 440 LB	
NET WEIGHT 2.3 KG	DIVISION 0.1 KG / 0.2 LB	UNIT KG / LB / ST:LB	POWER SUPPLY AA x 4
FUNCTION LOW POWER & OVERLOAD INDICATION AUTO POWER OFF IN 15 SECONDS			

FEATURES

- NO need to download APP or pair devices
- Flexible and safe data management
- Automatically upload data
- Remote care anytime, anywhere
- Durable tempered glass surface



CELLULAR INFRARED THERMOMETER

4G
NETWORK



JPD-FR409

DISPLAY	MEMORY STORAGE	TEMPERATURE UNITS
HD LCD DISPLAY	20 SETS OF MEASUREMENT DATA	°C/°F ADJUSTABLE
BATTERY	AUTOMATIC SHUTDOWN FUNCTION	
AAA X 2	10S ± 1S	
MEASURING RANGE & ACCURACY		
FOREHEAD MODE	22.0°C - 43.0°C (71.6°F-109.4°F)	
EAR MODE	34.0°C-43.0°C (93.2°F-109.4°F)	

FEATURES

- NO need to download APP or pair devices
- Flexible and safe temperature management
- Automatically upload data
- Remote care anytime,anywhere

