



Official Website



JUMPER Health

Shenzhen Jumper Medical Equipment Co., Ltd.

Address: D Building, No. 71, Xintian Road, Fuyong Street, Baoan, Shenzhen, Guangdong, China

Tel: +86 755-26696279 **E-mail:** info@jumper-medical.com **Web:** <https://jumpermed.com>

PRODUCT BROCHURE

SINCE
1998




JUMPER

Healthcare within mind


<https://jumpermed.com>

POCKET FETAL DOPPLER



HEAR THE BABY'S HEARTBEAT.

Pocket Fetal Doppler is a hand-held fetal heart rate detecting device. It displays fetal heart rate on a digital LCD screen /TFT screen with easy and convenient operation, which can be used in hospitals, clinics and home for daily self-check.




SHB10

DISPLAY LED	ULTRASOUND FREQUENCY 3.0 MHz	BATTERY AAA X 5
-----------------------	--	---------------------------

NET WEIGHT 156.5g (WITHOUT BATTERY)	DIMENSIONS(mm) 154 x 69 x 32	FHR MEASUREMENT RANGE 50-240 BPM
--	--	--

PERFORMANCE BATTERY STATUS FHR VALUE VOLUME	MEASURING MODE REAL-TIME	FUNCTION ADJUSTABLE VOLUME
---	------------------------------------	--------------------------------------




SHB20

DISPLAY LED	ULTRASOUND FREQUENCY 3.0 MHz	BATTERY AAA X 5
-----------------------	--	---------------------------

NET WEIGHT 156.5g (WITHOUT BATTERY)	DIMENSIONS(mm) 154 x 69 x 32	FHR MEASUREMENT RANGE 50-240 BPM
--	--	--

PERFORMANCE BATTERY STATUS FHR VALUE VOLUME	MEASURING MODE REAL-TIME	FUNCTION ADJUSTABLE VOLUME
---	------------------------------------	--------------------------------------




JPD-100F

DISPLAY TFT	ULTRASOUND FREQUENCY 2.5 MHz	BATTERY 3.7 V Lithium (RECHARGEABLE)
-----------------------	--	---

NET WEIGHT 250g	DIMENSIONS(mm) 139 x 82 x 28	FHR MEASUREMENT RANGE 50-210 ± 2 BPM
---------------------------	--	--

PERFORMANCE BATTERY STATUS WAVE FORM FHR VALUE VOLUME	MEASURING MODE REAL-TIME	FUNCTION ADJUSTABLE VOLUME
---	------------------------------------	--------------------------------------





SHA10

DISPLAY TFT	ULTRASOUND FREQUENCY 3.0 MHz	BATTERY AA x 2
NET WEIGHT 130g	DIMENSIONS(mm) 142.6 x 40.5 x 42.5	FHR MEASUREMENT RANGE 50-210 ± 2 BPM
PERFORMANCE BATTERY STATUS WAVE FORM FHR VALUE VOLUME	MEASURING MODE REAL-TIME	FUNCTION ADJUSTABLE VOLUME BACKLIGHT
		
FHR	Clinically Tested	From 9th Week
		
Battery Included	TFT Screen	1 min Auto Power Off



JPD-100E2
(AAA X 3)

JPD-100E

DISPLAY TFT	ULTRASOUND FREQUENCY 2.5 MHz	BATTERY 3.7 V Lithium (RECHARGEABLE)
NET WEIGHT 250g (WITHOUT BATTERY)	DIMENSIONS(mm) 147 x 101.9 x 31.6	FHR MEASUREMENT RANGE 50-210 ± 2 BPM
PERFORMANCE BATTERY STATUS WAVE FORM FHR VALUE VOLUME	MEASURING MODE REAL-TIME	FUNCTION ADJUSTABLE VOLUME
		
FHR	Clinically Tested	From 16th Week
		
Battery Included	TFT Screen	1 min Auto Power Off



SHA20

DISPLAY TFT	ULTRASOUND FREQUENCY 3.0 MHz	BATTERY AA x 2
NET WEIGHT 130g	DIMENSIONS(mm) 142.6 x 40.5 x 42.5	FHR MEASUREMENT RANGE 50-210 ± 2 BPM
PERFORMANCE BATTERY STATUS WAVE FORM FHR VALUE VOLUME	MEASURING MODE REAL-TIME	FUNCTION ADJUSTABLE VOLUME BACKLIGHT
		
FHR	Clinically Tested	From 9th Week
		
Battery Included	TFT Screen	1 min Auto Power Off



JPD-100A2
(AA X 3)

JPD-100A

DISPLAY HD LCD	ULTRASOUND FREQUENCY 3.0 MHz	BATTERY 7.2V NI-MH (RECHARGEABLE)
NET WEIGHT 240g (WITHOUT BATTERY)	DIMENSIONS(mm) 135 x 83 x 33	FHR MEASUREMENT RANGE 50-210 ± 2 BPM
PERFORMANCE BATTERY STATUS FHR VALUE VOLUME	MEASURING MODE MANUAL COUNTING REAL-TIME AVERAGE	FUNCTION ADJUSTABLE VOLUME BACKLIGHT
		
FHR	Clinically Tested	From 16th Week
		
Battery Included	LCD Screen	1 min Auto Power Off



JPD-100E1
(AAA X 3)

JPD-100E+

DISPLAY TFT	ULTRASOUND FREQUENCY 2.5 MHz	BATTERY 3.7 V Lithium (RECHARGEABLE)
NET WEIGHT 250g	DIMENSIONS(mm) 147 x 101.9 x 31.6	FHR MEASUREMENT RANGE 50-210 ± 2 BPM
PERFORMANCE BATTERY STATUS WAVE FORM FHR VALUE VOLUME	MEASURING MODE REAL-TIME	FUNCTION ADJUSTABLE VOLUME
		
FHR	Clinically Tested	From 16th Week
		
Battery Included	TFT Screen	1 min Auto Power Off



JPD-100S6+

DISPLAY HD LCD	ULTRASOUND FREQUENCY 3.0 MHz	BATTERY AA x 2
NET WEIGHT 250g	DIMENSIONS(mm) 130 x 100 x 36	FHR MEASUREMENT RANGE 50-210 ± 2 BPM
PERFORMANCE BATTERY STATUS FHR VALUE	MEASURING MODE MANUAL COUNTING REAL-TIME AVERAGE	FUNCTION ADJUSTABLE VOLUME BACKLIGHT
		
FHR	Clinically Tested	From 16th Week
		
Battery Included	LCD Screen	1 min Auto Power Off



JPD-100S6-2

DISPLAY	ULTRASOUND FREQUENCY	BATTERY
TFT	3.0 MHz	AA x 2

NET WEIGHT	DIMENSIONS(mm)	FHR MEASUREMENT RANGE
250g	130 x 100 x 36	50-210 ± 2 BPM

PERFORMANCE	MEASURING MODE	FUNCTION
BATTERY STATUS WAVE FORM FHR VALUE	MANUAL COUNTING REAL-TIME AVERAGE	ADJUSTABLE VOLUME BACKLIGHT



FHR



Clinically
Tested



From
16th Week



Battery
Included



TFT
Screen



1 min Auto
Power Off



JPD-100B3
(AA X 2)

JPD-100B4

DISPLAY	ULTRASOUND FREQUENCY	BATTERY
TFT	2.5 MHz	9.6V NI-MH (RECHARGEABLE)

NET WEIGHT	DIMENSIONS(mm)	FHR MEASUREMENT RANGE
240g (WITHOUT BATTERY)	129 x 83 x 34	50-210 ± 2 BPM

PERFORMANCE	MEASURING MODE	FUNCTION
BATTERY STATUS WAVE FORM FHR VALUE VOLUME	REAL-TIME	ADJUSTABLE VOLUME BACKLIGHT



FHR



Clinically
Tested



From
16th Week



Battery
Included



TFT
Screen



1 min Auto
Power Off



JPD-100S6

DISPLAY	ULTRASOUND FREQUENCY	BATTERY
HD LCD	3.0 MHz	AA x 2

NET WEIGHT	DIMENSIONS(mm)	FHR MEASUREMENT RANGE
250g	130 x 100 x 36	50-210 ± 2 BPM

PERFORMANCE	MEASURING MODE	FUNCTION
BATTERY STATUS FHR VALUE	MANUAL COUNTING REAL-TIME AVERAGE	ADJUSTABLE VOLUME BACKLIGHT



FHR



Clinically
Tested



From
16th Week



Battery
Included



LCD
Screen



1 min Auto
Power Off



JPD-100B+

DISPLAY	ULTRASOUND FREQUENCY	BATTERY
HD LCD	2.5 MHz	7.4V Lithium (RECHARGEABLE)

NET WEIGHT	DIMENSIONS(mm)	FHR MEASUREMENT RANGE
250g (WITHOUT BATTERY)	195 x 190 x 52	50-210 ± 2 BPM

PERFORMANCE	MEASURING MODE	FUNCTION
BATTERY STATUS FHR VALUE	MANUAL COUNTING REAL-TIME AVERAGE	ADJUSTABLE VOLUME BACKLIGHT



FHR



Clinically
Tested



From
16th Week



Battery
Included



LCD
Screen



1 min Auto
Power Off



JPD-100B2
(AA X 2)

JPD-100B

DISPLAY	ULTRASOUND FREQUENCY	BATTERY
HD LCD	2.5 MHz	9.6V NI-MH (RECHARGEABLE)

NET WEIGHT	DIMENSIONS(mm)	FHR MEASUREMENT RANGE
240g (WITHOUT BATTERY)	129 x 83 x 34	50-210 ± 2 BPM

PERFORMANCE	MEASURING MODE	FUNCTION
BATTERY STATUS FHR VALUE	REAL-TIME	ADJUSTABLE VOLUME BACKLIGHT



FHR



Clinically
Tested



From
16th Week



Battery
Included



LCD
Screen



Adjustable
Volume

Pocket Fetal Doppler			
Parameters			
 Product Model	JPD-100B/JPD-100B2	JPD-100B3/JPD-100B4	
 Ultrasound Frequency	2.5 MHz		
 Display Screen	LCD	TFT	
 Bluetooth	—		
 Volume Display	—	✓	
 Wave Form	—	✓	
 Measuring Mode	① Real-time		
 Battery Type	9.6V NI-MH / AA x 2 (RECHARGEABLE)		
 Dimensions (mm)	129 x 83 x 34		
 Weight(g) (Include Batteries)	240 (WITHOUT BATTERY)		
 Features	① FHR Measurement Range: 50-210 ± 2 bpm ② Volume Adjustable ③ Battery Status ④ 1 min auto power off		